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
WELCOME GUIDE

MagniSonic | OPC-ET500N |

*Please read the instructions carefully before using the product and keep the user manual in a safe place for future reference.


PACKAGE CONTENTS & PRODUCT OVERVIEW

BASIC PRODUCT SPECIFICATIONS

Cover Material: PC + ABS
Battery Capacity: 500mAh
Normal Usage Time: ≥ 60 days
Charging Time: ≤ 6 hours
Rated Working Voltage: DC 3.7V
Rated Working Power: 2W
*oraimo Lab test results

INSTALLATION INSTRUCTIONS

Step 1: Attach the Brush Head
Align the brush head with the handle and press it down firmly onto the metal shaft.

Step 2: Ensure a Secure Fit
Gently pull up and push down the brush head to check if it is firmly attached.

Note:
Ensure the notch on the neck aligns with the inner notches of the toothbrush head to lock it in place.

USAGE INSTRUCTIONS

Oral Cavity

1. Wet the bristles and apply a suitable amount of toothpaste. To prevent splattering, position the toothbrush on your teeth before turning it on.
2. Angle the bristles slightly towards the gum line and position them against your teeth, hold the brush at a 45° angle to your teeth and gums.
3. Press the power switch and gently glide the brush head across your teeth using small, vertical motions. For each section (as illustrated), brush the front, back, and chewing surfaces for 10 seconds each.
4. After brushing, remove the toothbrush from your mouth, rinse it thoroughly, and store it in a well-ventilated, dry location.

MAIN SPECIFICATIONS

Brushing Mode

Gear Mode	Frequency/Speed	Description	Product Status
Gear 1	32000 times/minute	Strong gear	The first white indicator light is on
Gear 2	32000 times/minute	Mid-range	The second white indicator light is on
Gear 3	32000 times/minute	Weak gear	The third white indicator light is on
Gear 4	32000 times/minute	Massage mode	The fourth white indicator light is on
Gear 5	36000-47000 times/minute	Variable speed	The dynamic frequency will automatically, and the fifth white indicator light is on

Battery Indicator

Battery Level	Indicator Light After Shutdown	Product Status
Less than 10%	No indicator light	Product not working
10% - 20%	One gear light	Product in low power mode
20% - 40%	Two gear lights	Product in medium power mode
40% - 60%	Three gear lights	Product in medium-high power mode
60% - 80%	Four gear lights	Product in high power mode
Greater than 80%	Five gear lights	Product fully charged

Charge Indicator

Battery Level	Indicator Light Status
Less than 10%	Charging indicator light not on
10% - 90%	Gear lights gradually light up from bottom to top, indicating a cycle
90% - 100%	Charging indicator light not on

TROUBLESHOOTING

Trouble	Solutions
The product does not work.	Check if the battery has power.
The vibrational frequency has decreased.	The battery may be low. Please charge it in time.
The brush head has less cleaning effect.	Replace the toothbrush head at least once every three months.

IMPORTANT SAFETY INSTRUCTIONS

1. Do not open, repair, or modify the product without authorization.
2. Ensure toddlers or individuals unable to operate the device are supervised by an adult.
3. Avoid submerging the product in water for extended periods.
4. Remove the battery if the product will not be used for a long time.
5. Do not exceed the product's intended functional use.
6. This product is designed for household use to clean teeth, gums, and tongue only. Do not use it for other purposes.
7. Use the vibrational strength provided by the toothbrush; do not apply excessive pressure on the teeth or gums.
8. Keep the product away from high temperatures and direct sunlight.
9. Avoid striking or applying excessive force to the product.

WARRANTY CERTIFICATE

Name :	
Address :	
Tel. No. :	
Date of purchase : (dd/mm/yy)	
Store name :	
Warranty period :	

WARRANTY TERMS

1) The warranty period starts from the date of delivery of the product and is 2 years. (This period may not be less than 2 years.)
2) The product, including all its parts, is covered under warranty.
3) If the product is found to be defective, the consumer may exercise one of the rights specified in Article 11 of Law No. 6502 on the Protection of Consumers, namely:
 a- Withdrawal from the contract,
 b- Requesting a price reduction,
 c- Requesting free repair,
 d-Requesting replacement with a non-defective equivalent.
4) If the consumer chooses the right of free repair, the seller is obliged to repair or have the product repaired free of charge, without demanding any fee under any name such as labor cost, replaced part cost, or otherwise. The consumer may also exercise this right against the manufacturer or importer. The seller, manufacturer, and importer are jointly and severally liable.
5) If the consumer uses the right of free repair and:
 - The product breaks down again within the warranty period,
 - The maximum repair period is exceeded,
 - It is determined by an authorized service station, seller, manufacturer, or importer (by report) that repair is not possible;
 the consumer may request a refund, a price reduction proportional to the defect, or -replacement with a non-defective equivalent, if possible. The seller may not refuse this request. If the request is not fulfilled, the seller, manufacturer, and importer are jointly and severally liable.
6) The repair period shall not exceed 20 business days. This period starts from the date the fault is reported to the authorized service station or seller within the warranty period, or from the date the product is delivered to the authorized service station if outside the warranty period.
 If the defect cannot be repaired within 10 business days, the manufacturer or importer must provide the consumer with another product of similar features for use until the repair is completed.
7) Failures resulting from use contrary to the instructions in the user manual are outside the scope of the warranty.
8) In disputes arising from the exercise of warranty rights, the consumer may apply to the Consumer Arbitration Committee or the Consumer Court in the place of residence or where the consumer transaction was carried out.
9) If this Warranty Certificate is not provided by the seller, the consumer may apply to the General Directorate of Consumer Protection and Market Surveillance of the Ministry of Trade.

Technical Service
Ouno Servis Tuzla Şube
Address: Tepeören Demokrasi Caddesi No:137 ,34959 Tuzla
Sadece kargo gönderimlerine hizmet vermektedir.
Ouno Servis Genel Merkez
Address: Atilla İlhan Cd. No:20, 34750 Ataşehir/İstanbul
Sadece ön bürodan hizmet verilmektedir.
Website: https://ounoservice.com/
For authorized service centers, please visit https://ounoservice.com/, as required by the distributor.

Importer Company:
Company Name: Dora Communication and Trade Inc.
Address: Küçükbakkalköy Mah., Atilla İlhan Cd. No:22, 34750 Ataşehir, Istanbul
Telephone: (0216) 469 40 60
Authorized Signature:
Company Stamp:

Manufactured for:
Company Name: oraimo TECHNOLOGY LIMITED
Address: FLAT N 16/F BLOCK B UNIVERSAL INDUSTRIAL CENTRE 19-25 SHAN MEI STREET FOTAN NT HONGKONG
Email: hello.ab@oraimo.com

Manufactured by:
Company Name: oraimo TECHNOLOGY LIMITED
Address: FLAT N 16/F BLOCK B UNIVERSAL INDUSTRIAL CENTRE 19-25 SHAN MEI STREET FOTAN NT HONGKONG
Email: hello.ab@oraimo.com

Transport and Shipping Instructions

- Before moving the device, ensure all connections are disconnected and cables removed.
- Use the original box and packaging materials for repacking.
- Avoid shock, impact, heat, humidity, and dust during use and transport.

Consumer Maintenance and Repairs
No consumer serviceable or repairable parts. For maintenance or repair, contact an authorized service point listed in the Technical Service List.

Efficient Energy Use
Turn off data connections (Bluetooth, Wi Fi) when not in use.
Lower screen brightness and set shorter screen timeout.
Disable unnecessary sounds

Recycling

This product complies with the AEEE Regulation.

Proper disposal of the product.
This symbol indicates that the product should not be disposed of together with other household waste across the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal and to promote the sustainable reuse of material resources, please recycle responsibly. To return a used product, please use return and collection systems or contact the retailer from whom the product was purchased. They can take this product for environmentally safe recycling.

SCAN↓DOWNLOAD

oraimo store
Easily Smart, Simply Best

Authentic Products | Fast Delivery | Best Price | Prompt Service

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