

Your Manufacturer and Master Distributor for Physiotherapy, Sports & Fitness Products

Version 2021/10

We are proud to present to you our new catalogue 2021-2022!

Since 1996, our company has been known as a worldwide manufacturer and master distributor of a broad range of physiotherapy, rehabilitation and sports & fitness products. Our products reach the end-users through a global network of dealers and distributors.

Our own branded MoVeS and Mambo Max products, as well as other well-known brands in our market, are all presented in this catalogue.

On page 7, you can find our new products: our MoVeS Kinesiology Tape, Mambo Max Balance Pad Pro, and our updated MoVeS PearlPack Sport Pack with strap.

All our medical devices, marked with the MD symbol, comply with all the requirements of the EU's Medical Devices Regulation (MDR), which came into force May 26, 2021.

For all your inquiries for Private Label products, please contact us. We will be happy to assist you working out your private label projects.

Unfortunately, Covid-19 is still around, which has made 2021 a challenging year as we are currently experiencing the economic impact of this pandemic. Shortages of raw materials & containers as well as available vessel spaces cause high product price and transport price increases, as well as many and long delays on deliveries. We have adapted our systems to this present-day reality to minimise stock shortages.

This year is also special in the company's history; I herewith take the opportunity to introduce you to the new owners. It is with great confidence that I hand over the company to Evelyn and Tom, who will continue in the same spirit and with enormous enthusiasm. I wish them lots of success!

Finally, I would of course like to thank YOU for your confidence in the company and its products! Stay safe, stay healthy!



01 - Resistive 8









MVS In Motion | Westdijk 150, 2830 Tisselt (Willebroek) | Belgium - Europe [T] +32 3 866 40 76 | www.mvs-in-motion.com | info@mvs-in-motion.com

<u>CONTENT</u>



///







TABLE OF CONTENT ///

OUR OWN BRANDS ONLINE	6 6	3. UPPER & LO
LEGEND	6	Fit Stretch
NEW PRODUCTS	7	Saehan Cycle
THE ENVIRONMENT	7	Wheeleo
		MoVeS Oxycy
1. RESISTIVE	11	MoVeS Oxycy
		MoVeS Oxycy
MoVeS Band	11	MoVeS Oxycy
MoVeS Latex-Free Band	13	Curamotion
MoVeS Tube	14	MoVeS Pedal
MoVeS Wide Loop	15	MoVeS Fold
MoVeS Loop	15	
MoVeS F!T Band	17	MoVeS Stan
MoVeS F!T Tube	18	MoVeS Stan
MoVeS F!T Loop	19	MoVeS Wood
MoVeS F!T Superloop	19	MoVeS Woo
MoVeS Bar	20	MoVeS Woo
MoVeS Ring Tube	20	Plastic Incline
Multi-Grip Exerciser	21	
MoVeS Ankleciser	21	
MoVeS Shoulder Rope Pulley - Metal Bracket	22	4. WELLNESS
MoVeS Shoulder Tube Pulley	22	
MoVeS Knee Rope Pulley	22	Kinesiology Ta
MoVeS Shoulder Rope Pulley - Door Anchor	23	Yoga
Notch Handles	23	Stepping Ston
Value-Pack	23	Triggerpoint R
	_	Spiky Foam
Wooden Handles	23	Hollow Foar
Hand Support	23	Ergonomic I
MoVeS Accessories	24	3-in-1 Foam
		Massage Sti
2. HAND	25	Sitting Cushio
		Pilates
MoVeS Flex-Ion	25	Soft-Over-B
Digi-Flex	26	Foam Roller
Digi-Extend	26	Ring
Prohands Series	26	Massage Tool
Handmaster Plus	27	Massage Ball,
Constant Force X-Tend	27	
Varigrip Therapy	28	
Varigrip Sport	28	5. FITNESS
Varigrip Uno	28	61
Power-Web Combo	29	Gym Mat
Power-Web Flex-Grip	29	Ball Accessori
Power-Web Senior	30	Aerobic Step
Power-Web Junior	30	Fabric Elastic
MoVeS Squeeze Balls & Eggs	31	Dynadome
MoVeS Wrist Exerciser	31	Balance Board
Adjustable Hand Grip	32	Balance Board
Foam Hand Grip	32	Wooden Bala
MoVeS Comfort Putty	33	AB Gym Ball
Theraflex Putty	33	AB Peanut Ba
		Body Twister
		Ekilibre

3.	UPPER & LOWER BODY	34
	Fit Stretch	34
	Saehan Cycle	34
	Wheeleo	35
	MoVeS Oxycycle 1 - Active Magnetic	36
	MoVeS Oxycycle 2 - Passive	36
	MoVeS Oxycycle 3 - Active & Passive	37
	MoVeS Oxycycle Pedals	37
	Curamotion	38
	MoVeS Pedal Exerciser	39
	MoVeS Foldable Pedal Exerciser	39
	MoVeS Standard Pedal Exerciser Premium Chrome	39
	MoVeS Standard Pedal Exerciser K/D	39
	MoVeS Wooden Incline Board	40
	MoVeS Wooden Incline Board - Ladder	40
	MoVeS Wooden Incline Board - Fix	40
	Plastic Incline Board	40
		40
4.	WELLNESS	41
	Kinesiology Tape	41
	Yoga	42
	Stepping Stones	42
	Triggerpoint Rollers	43
	Spiky Foam Roller	43
	Hollow Foam Roller	43
	Ergonomic Foam Roller	43
	3-in-1 Foam roller	43
	Massage Stick	43
	Sitting Cushions	44
	Pilates	45
	Soft-Over-Ball	45
	Foam Roller	45
	Ring	45
	Massage Tools	46
	Massage Ball, Roll & Ring	46
5.	FITNESS	47
	Gym Mat	47
	Ball Accessories	48
	Aerobic Step	48
	Fabric Elastic Band	49
	Dynadome	49
	Balance Board Freeman	50
	Balance Board	50
	Wooden Balance Boards	50
	AB Gym Ball	51
	AB Peanut Ball	51
	Body Twister	52
	Ekilibre	52
	Wooden Work-Out Station	52
	Balance Pad Pro	53
	Balance Pad	53
	Balance Trainer Swim Noodle	53 54
	Swim NUULIE	54

FAIRS	54
6. STRENGTH	55
Coated Dumbbell	55
SoftMeds	55
Wrist & Ankle Weight	56
Thumb Lock Wrist Weight	56
Medicine Ball	57
Kettlebell	57
7. THERMAL	58
Para Pro - Paraffin Bath	58
PearlPack	59
PearlPack Knee/Elbow	59
PearlPack Eye Mask	59
PearlPack Back Wrap	59
PearlPack Sports Pack + strap	59
Roll-On Pack	60
Torex Pro Ankle Boot	60
Hot/Cold Pack	61
8. EVALUATION	62
Infrared Thermometer	62
Ear & Forehead Thermometer	62
Shotgun Thermometer	62
Blood Pressure Monitor	63
Fingertip Pulse Oximeter	63
Adult Oximeter	63
Junior Oximeter	63
Hydraulic 3-piece hand evaluation kit	64
Hydraulic 7-piece hand evaluation kit	64



Hydraulic hand dynamometer	65
Digital hand dynamometer with Data Transfer Software-	
DHD-1	65
Digital hand dynamometer with G-STAR™ Software -	
DHD-3	65
Smedley Hand Dynamometer	66
Push-Pull Dynamometer	66
Body Fat Scale	66
Pinch Gauge	67
Squeeze Dynamometer	67
Back-Leg-Chest Dynamometer	67
Fat Caliper	68
Measurement Tape	68
Tuning Fork	69
Inclinometer	69
Scoliometer	69
Goniometers	70
Baseline	70
Lafayette	70
Goniometer Kit	70
Plastic Goniometer	71
Stainless Steel Goniometer	71
Finger Goniometer	71
Monofilament	72
Discriminator	72
Neurological Reflex Instruments	73
Sit & Reach Unit	73
TENS device	74
PRIVATE LABEL	75
AFTER-SALES	76
INDEX	77
TERMS & CONDITIONS	78

OUR OWN BRANDS ///





MoVeS products are aimed towards physiotherapy and rehabilitation, predominantly resistive exercise. Resistive exercise products have various levels of resistance, which also makes them perfect for fitness. MoVeS products are generally MD certified.




Our Mambo Max products are the go-to-products for wellness, fitness and strength training at an affordable price tag. Combine various products for a complete work-out, without leaving your home.

ONLINE



///

Kinesiology Tape P41 MoVeS NEW





We, at MVS In Motion, continue our strive to reduce our carbon footprint on the planet. Because after all, there is no planet B!

By making exercises available on our websites, not only do we reduce our ecological footprint by cutting down on the number of unused prints, it also allows us to keep costs low and offer you better value-for-money.



www.moves-you.com www.moves-band.com

The virtual home of our MoVeS and MoVeS F!T product line. This B2C platform connects our distributors directly to the end users, with product information, a shop locator and numerous free downloadable exercises. The website has been upgraded and is mobile friendly.

www.mambomax.com

Completely revitalised! This is the virtual home of our Mambo Max product line. This B2C platform connects our products directly to the end users, with product information and numerous free downloadable exercises.

www.mvs-in-motion.com

MVS In Motion's B2B site, with all available products and digital catalogue. Ask for your login to download all high resolution pictures, product information, logos, videos and many more marketing tools.

Social Media in 🖸 f 💿

Here you will find our latests and useful information. Like or follow our pages for a higher brand awareness to create more demand. Our pages are continuously monitored to answer your questions as soon as possible. We are active on LinkedIn and YouTube. Soob, you will also

be able to find us on Facebook and Instagram




In order to keep all the information in this catalogue as clear and easy accessible as possible, we use a system of icons that are shown next to the product where applicable.

MD

	×
	the second
m	





instructions.





Packed in a polybag with full colour header.



Plastic free product and packaging

Contains multilingual

Packed in a blister pack with full colour insert.

Packed in a shrinkfoil wrap.	Ø	Product Diameter
Comes with a pump.	\$	Product Height
Packaging contains printed exercises.	\leftrightarrow	Product Width
Downloadable exercises available on our websites.	R	Maximum User We
Medical Device, in compliance with MDR 2017/745	CE	European Conformi (Health, safety & environment)

۶	
\$	Product Height
↔	Product Width
ß	Maximum User Weight
((European Conformity

NEW PRODUCTS

///

THE ENVIRONMENT ///



We have also been making continuous efforts over the past years to replace our plastic packaging with the more eco-friendly and recyclable cardboard packaging. And we will continue these efforts in the future!

MoVeS



LEADER IN TECHNOLOGY

What is Snap-Stop[®]?

Snap-Stop^{*} makes exercising with latex bands much safer. This unique feature will avoid sudden snapping, even if there's a cut or puncture. When working out with a resistive exercise band, the band has to endure extremely high forces. A small cut can have a big consequence...

MoVeS and MoVeS F!T Bands are unique in their field. Thanks to our extensive research and development we have achieved an innovative production process: when used properly, the latex band will not suddenly snap, even when exposed to high forces. This does not apply to a very old band, or band that has been exposed outdoors for some time, as well as an abused band. Using a fresh and clean sheet of band is recommended.

Snap-Stop[®] for safer resistive exercising.

1000


One-Stop-Shop for all your bands and tubes

MVS In Motion is your one-stop-shop for all resistive exercise products. With MoVeS we offer both original colour and F!T colour sequence. In addition, we offer a wide range of Latex-Free resistive exercise bands in original colours. Our tubes and loops are also available in the original and F!T colour sequence.

The latex bands, tubes and (super)loops are made of high quality latex from huge rubber tree plantations in Malaysia. All products have an extraordinary elasticity and linearity. The Snap-Stop® feature on our latex bands, dipped production process for our tubes and the multi-layered structure of our (super)loops makes them your ideal partner for your physiotherapy trainings or intensive work-outs. The different resistances allow you to use MoVeS for rehabilitation as well as for fitness.

MoVeS Latex-Free is made from TPE, to the same high quality standards as the regular MoVeS bands. Never before were the characteristics from Latex-Free bands so close to those from their latex counterpart.

MoVeS Band: 8 resistances (Tan - Gold) & 5 lengths (1,5m - 2,5m - 5,5m - 22,5m - 45,5m) **MoVeS F!T Band:** 6 resistances (Peach - Grey) & 4 lengths (1,5m - 2,5m - 5,5m - 45,5m) **MoVeS Latex-Free Band:** 5 resistances (Yellow - Black) & 3 lengths (1,5m - 5,5m - 45,5m)






- Well-known colours used by physiotherapists all over the world
- Packaging with two tone layout brings more variety
- MoVeS brings premium quality to an affordable price range.












Resistance - MoVeS Band vs. Market Leader





MoVeS Bands can be stretched up to 600% elongation! MoVeS Bands are therefore your perfect partner for a rehabilitation programme or a full body work-out. Even after an effective training with amazing results, MoVeS Bands allow you to use the same band again. When you prestretch your band, you will be able to achieve higher forces without having to worry about the 'elastic stop', the point where the band feels like a regular string. Most competitor's bands are not to be stretched more than 3 times their resting length, while a MoVeS Band achieves easily the double, maintaining its amazing elasticity.

These specifications apply to all MoVeS resistances. The high-quality material enables you to have a much wider array of resistances and possibilities compared to the competitors. At the start of your exercise, the band already gives some resistance, which allows you to achieve much better results than when the band is at its resting length.

Note:

- 50% elongation = 1,5 x lay-flat length
- 100% elongation = 2 x lay-flat length
- 150% elongation = 2,5 x lay-flat length

cb po m ex ex

MoVeS Bands are premium quality, low powder (in fact, less than 2mg/dm², so we're officially powder-free) latex resistive exercise bands. The unique Snap-Stop* fabrication process makes this product extremely tear resistant. Using the MoVeS Bands for resistive exercises provides positive as well as negative forces on the muscles and improves the range of motion, the strength and the cooperation of the muscle groups. The 8 different resistance levels are determined by the thickness of the band. MoVes bands are available in 5 lengths and come with a full colour manual. The band has a width of +/- 127mm.

Dispenser 25 bands of 1,5m in a nicely designed OTC dispenser box. Each band comes in a separate polybag.

- 1,5m 10 pre-cut bands of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header and barcode.
- 2,5m 10 pre-cut pieces of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header and barcode.
- 5,5m One small roll of our premium band.
- 22,5m One medium roll of our premium band.
- 45,5m One large roll of our premium band.



				MOVES	BAND		
COLOUR	RESISTANCE	1,5m Dispenser 25pcs	1,5m 10-pack	2,5m 10-pack	5,5m Box Small Roll	22,5m Box Medium Roll	45,5m Box Large Roll
Tan	Extra Light	n/a	n/a	n/a	01-100501	n/a	01-104501
Yellow	Light	01-100102	01-101502	01-100202	01-100502	01-102202	01-104502
Red	Medium	01-100103	01-101503	01-100203	01-100503	01-102203	01-104503
Green	Heavy	01-100104	01-101504	01-100204	01-100504	01-102204	01-104504
Blue	Extra Heavy	01-100105	01-101505	01-100205	01-100505	01-102205	01-104505
Black	Special Heavy	01-100106	01-101506	01-100206	01-100506	01-102206	01-104506
Silver	Super Heavy	n/a	01-101507	01-100207	01-100507	01-102207	n/a
Gold	Ultra Heavy	n/a	n/a	n/a	01-100508	01-102208	n/a
	5 resistances w to black)	n/a	01-101511 (5-pack)	01-100211 (5-pack)	n/a	n/a	ap-st
						Sr	Jar

MoVeS Band $MD \in$



	TOTAL KG / ELONGATION %					
COLOUR -	50%	100%	150%	200%	250%	
Tan	0,8	1,1	1,6	1,9	2,2	
Yellow	1,0	1,3	1,5	1,9	2,3	
Red	1,5	2,0	2,5	3,2	3,7	
Green	1,6	2,5	3,0	3,5	4,1	
Blue	2,2	3,0	3,6	4,4	5,2	
Black	2,5	3,4	4,0	4,9	6,1	
Silver	4,1	5,8	7,2	7,9	9,5	
Gold	5,8	7,9	10,0	12,0	14,0	


02-HAND

11

Resistance - MoVeS Latex-Free Band vs. Market Leader –

Research and development is high on our list for MoVeS products. As an example, these continuous efforts have brought us the Snap-Stop feature for the regular MoVeS Band.

Now, we have improved our Latex-Free Bands. The MoVeS Latex-Free Band has been optimised by using an innovative mix of raw material. Following customer demand, the grip has been improved. Resistance levels are similar to the 'old' latex-free band.



We are confident that MoVeS Latex-Free Bands are one of the best latex-free resistive exercise bands on the market. The elasticity and linearity are simply amazing. Finally people who don't support latex can exercise with the same high quality specifications. MoVeS Latex-Free bands start from high quality thermoplastic elastomer (TPE), so it's non-toxic and completely PVC and phthalate free.

The linearity from how the strength per elongation evolves is one of the most important specifications as it offers the possibility to apply a very extended range of exercises without risking an injury due to strange unexpected product behaviour. Latex-free alternatives will reach their 'elastic stop' more quickly; this is the point where the band feels like a regular string.

We notice a growing demand and trend for latex free products. For instance, more and more hospitals and centers aim to be completely latex free to avoid possible allergic reactions for their patients or staff. After extensive research and development, we were able to launch our MoVeS Latex-Free Bands. A product which we are extremely proud of.



Thanks to our technological leadership, we achieve a Latex-Free exercise band with characteristics that were never before so close to its latex counterpart. With an exceptional 'resistance vs. elongation' linearity, we are confident our MoVeS Latex-Free Band is the best product to fulfil your needs.

The 5 different resistance levels are determined by the thickness of the band. These bands are available in 3 lengths and come with a full colour manual. The band has a width of +/- 100mm.

- 1,5m 10 pre-cut bands of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header.

		TOTA	L KG / EL	ONGATI	ON %	
COLOUR	50%	100%	150%	200%	250%	300%
Yellow	1,3	1,7	1,9	2,3	2,8	3,2
Red	1,4	1,9	2,2	2,6	3,0	3,5
Green	2,0	2,5	3,0	3,4	4,0	4,5
Blue	2,3	3,0	3,4	3,9	4,7	5,5
Black	2,4	3,4	3,9	4,3	5,1	6,1



			MOV	ES LATEX-FREE	BAND
	COLOUR	RESISTANCE	1,5m	5,5m Box	45,5n
_			10-pack	Small Roll	Large
	Yellow	Light	01-121502	01-120502	01-12
	Red	Medium	01-121503	01-120503	01-12
	Green	Heavy	01-121504	01-120504	01-12
	Blue	Extra Heavy	01-121505	01-120505	01-12
	Black	Special Heavy	01-121506	01-120506	01-12
	Set of 5	5 resistances	01-121511	n/2	n
_	(yello	w to black)	(5-pack)	n/a	n/





02-HA









The MoVeS Tube is a high quality tube created by an intensive dipping process. The dipping process guarantees a multi-layered tube, which ensures a strong and safe training product. Extruded tubes are very susceptible to damage and can easily snap, this is not the case for our multi-layered dipped tubes. In general, dipped tubes are even more durable than bands, especially when used on a coarse surface.

MoVeS Tubes are available as a ready to use product - 120cm length with handles straight out of a polybag - and in rolls - 7,5m or 30m packed in a nicely designed box.

The table below shows the exceptional strength to elongation linearity from our tubes.

		TOTA	l KG / El	ONGATI	ON %	
COLOUR -	50%	100%	150%	200%	250%	300%
Yellow	1,0	1,5	1,9	2,3	2,9	3,6
Red	1,4	2,0	2,5	2,8	3,4	4,0
Green	2,4	3,4	4,0	4,9	5,5	6,5
Blue	2,9	3,9	4,7	5,6	5,9	7,4
Black	3,4	4,4	5,4	6,2	7,1	8,1
Silver	4,0	5,4	6,5	7,7	8,1	9,5

Combine a roll with our Safety Handles for superior grip (p22).

MoVeS Tubes with Handles are very useful for physiotherapy and muscle training. These

high quality and very resistant tubes are created to train all muscles in your body and improve body flexibility. The lighter resistances are perfect for physiotherapy or a light work-out, while the highest resistances will even convince the more trained athlete.



cb ex



We offer wide loops (7,5cm) next to our regular loops (2,5cm) with the same resistances. Some users may prefer the wider loops for spreading the forces more evenly on their body parts, while others prefer the thicker narrow loops which won't roll-up while being stretched.

As your supplier for bands, tubes and loops, these wide loops fit perfectly in our range of resistive products to satisfy all users. Wide loops are also multi-layered and available in 6 resistances: from Extra Light (Tan) to Special Heavy (Black).

MoVeS Wide Loops:

High quality, cost effective and easy to use! Lay-flat-length: 30cm Width: 7,5cm One unit contains 10 individually packed pcs of the same colour/resistance

COLOUR	RESISTANCE	MOVES WIDE LOOP	UNIT
Tan	Extra Light	01-100071	10 pcs
Yellow	Light	01-100072	10 pcs
Red	Medium	01-100073	10 pcs
Green	Heavy	01-100074	10 pcs
Blue	Extra Heavy	01-100075	10 pcs
Black	Special Heavy	01-100076	10 pcs

MD (E MoVeS Loop **M**_oV_eS MoVeS Loops are extremely portable and unbelievably versatile products. They are very effective

instruments for stretching and strengthening all muscle groups. The MoVeS Loop Set really is a "Gym on-the-go".

go.

Gym on-the-90.

Our loops are also multi-layered which results in a more durable product. These thicker loops won't rollup during exercising, avoiding the uncomfortable feeling from thinner alternatives.

> They are available in 6 resistances: from Extra Light (Tan) to Special Heavy (Black). The loop set combines our most popular colours in a compact mesh bag. This affordable set is very useful to bring along wherever you

MoVeS Loops:

High quality, cost effective and easy to use! Lay-flat-length: 30cm Width: 2,5cm One polybag contains 10 loose pcs of the same colour/resistance

COLOUR	RESISTANCE	MOVES LOOP	UNIT	TOTAL KG / ELONGATION %						
Tan	Extra Light	01-100001	10 pcs	COLOUR	50%	100%	150%	200%	250%	300%
Yellow	Light	01-100002	10 pcs	Tan	2.2	3.2	3.8	4.6	5.4	6.4
Red	Medium	01-100003	10 pcs	Yellow	2.8	4,0	4,8	5,8	6,8	8,0
Green	Heavy	01-100004	10 pcs	Red	3,4	4,8	5,8	6,8	8,0	9,2
Blue	Extra Heavy	01-100005	10 pcs	Green	4,4	6,4	7,8	9,2	10,8	12,4
Black	Special Heavy	01-100006	10 pcs	Blue	5,0	7,2	9,0	10,8	12,6	14,4
Loop Set of 4 resistances		01-100011	Sot / A pcc	Black	6,1	8,9	11,0	13,0	15,0	17,3
(yellow to blue) in mesh bag		01-100011	Set / 4 pcs							

01-RESISTIVE

02

-HAND

03-UPPER & LOWER BODY

04-WELLNESS

07-THERMAI

14

MoVeS Wide Loop MD CE







PLASTIC
1 13333
FREE

Multi-Layered

01-RESISTIVE

02-HAND

03-UPPER & LOWER BODY

04-WELLNESS

Made in Malaysia 🦟



The production of MoVeS latex products start deep in the humid forests of Malaysia, where the best and strongest trees provide the natural latex rubber that will be processed from a milky fluid to your preferred resistive exercise band. When a tree is 5 - 6 years old, harvesting can begin: incisions are being made orthogonally in the bark of the tree, just deep enough to tap the latex vessels without harming the tree's growth. The fluid is collected in small buckets, a process known as rubber tapping.



Did you know that rubber trees combat 🖠 greenhouse gases and thereby slowing down global warming? Additionally, latex harvesting helps prevent deforestation.



natural **100%** biodegradable

"This is why we are confident when using the words 'premium quality' when referring to MoVeS."

Only the best quality latex is being selected for MoVeS Bands. Production starts under the constant supervision of a team of scientists, ensuring an immaculate product quality.

Latex sheets are produced, powdered, washed, dried and cut to the exact width. The amount of powder used is below the 2mg/dm² limit, which makes the MoVeS Bands officially powder-free. Each batch of latex sheeting that is produced, is subject to the meticulous inspection of our scientists by a series of laboratory experiments: colour-coding, dryness, tensile strength, tearing resistance, snapping resistance, ...

During the rolling process, a 100% visual inspection is being conducted on both sides of the band. Indeed, before leaving the factory, every MoVeS product has been completely analysed, checked and examined to ensure premium product quality.



MoVeS F!T, for athletes and home use.

cb po m ex ex

MoVeS F!T Bands are low powder (in fact, less than 2mg/dm², so we're officially powder-free) latex resistive exercise band with the same premium quality as the original MoVeS Bands. The difference is purely the alternative, more youthful colour sequence. The band comes in 6 different resistance levels which are determined by the thickness of the band, is available in 4 different lengths and comes with a full colour manual.

COLOUR -	T	OTAL KG	/ ELON	GATION	%	Compares to
COLOUR	50%	100%	150%	200%	250%	STANDARD
Peach	1,0	1,3	1,5	1,9	2,3	Yellow
Orange	1,5	2,0	2,5	3,2	3,7	Red
Lime Green	1,6	2,5	3,0	3,5	4,1	Green
Blueberry	2,2	3,0	3,6	4,4	5,2	Blue
Plum	2,5	3,4	4,0	4,9	6,1	Black
Grey	4,1	5,8	7,2	7,9	9,5	Silver

Our high quality 1,5m and 2,5m MoVeS F!T Bands are packed in retail packs. 10 pre-cut pieces of the same resistance or 5 pieces of different resistance that come in separate polybags with header for easy dispensing. Full-colour instruction manual is included. These 1,5m and 2,5m bands are easy to use and cutting is not required. The available colours are peach, orange, lime green, blueberry, plum and grey.





				MOVES F!T BAND)	
COLOUR	RESISTANCE	1,5m	2,5m	5,5m in Box	22,5m in Box	45,5m in Box
		10-pack	10-pack	Small Roll	Medium roll	Large Roll
Peach	Light	01-111502	01-110202	01-110502	n/a	01-114502
Orange	Medium	01-111503	01-110203	01-110503	n/a	01-114503
Lime Green	Heavy	01-111504	01-110204	01-110504	n/a	01-114504
Blueberry	Extra Heavy	01-111505	01-110205	01-110505	n/a	01-114505
Plum	Special Heavy	01-111506	01-110206	01-110506	n/a	01-114506
Grey	Super Heavy	01-111507	01-110207	01-110507	01-112207	P-Stop [®]
Set of 5 resistances (peach to plum)		01-111511	01-110211	n/a	n/	0-510
		(5-pack)	(5-pack)	11/a	Sna	il/a



17

MoVeS F!T Tube MD CE

The MoVeS FIT Tube is a high quality tube created by an intensive dipping process. The dipping process guarantees a multi-layered tube, which ensures a strong and safe training product. Extruded tubes are very susceptible to damage and can easily snap, this is not the case for our multi-layered dipped tubes. In general, dipped tubes are even more durable than bands, especially when used on a coarse surface.

MoVeS F!T Tubes are available as a ready to use product in 30m rolls, packed in a great dispenser box. Safety Handles offer an unlimited amount of exercises and superior grip.

MoVeS F!T Tubes are available in 6 different colours. The table below shows the resistance and to which standard colour it relates.

COLOUR -	٦	Compares to				
COLOOR	50%	100%	150%	200%	250%	STANDARD
Peach	1,0	1,5	1,9	2,3	2,9	Yellow
Orange	1,4	2,0	2,5	2,8	3,4	Red
Lime Green	2,4	3,4	4,0	4,9	5,5	Green
Blueberry	2,9	3,9		5,6	5,9	Blue
Plum	3,4	4,4	5,4	6,2	7,1	Black
Grey	4,0	5,4	6,5	7,7	8,1	Silver

Combine a roll with our
Safety Handles for superior grip (p22)

Compline a roll with our
Safety Handles for superior grip (p22).

Box of 30m contains maximal 3 separate pieces.



COLOUR	RESISTANCE	MOVES F!T TUBE
COLOUR	RESISTANCE	30m roll
Peach	Light	01-213002
Orange	Medium	01-213003
Lime Green	Heavy	01-213004
Blueberry	Extra Heavy	01-213005
Plum	Special Heavy	01-213006
Grey	Super Heavy	01-213007





MoVeS F!T Loops: High quality, cost effective and easy to use! Lay-flat-length: 30cm Width: 2,5cm

One Polybag contains 10 loose pcs



			and the second	
	COLOUR	RESISTANCE	MOVES F!T LOOP	UNIT
	Peach	Light	01-110002	10 pc
	Orange	Medium	01-110003	10 pc
	Lime Green	Heavy	01-110004	10 pc
- 2.	Blueberry	Extra Heavy	01-110005	10 pc
1		f 4 resistances perry) in mesh bag	01-110011	Set / 4

	·*.	The second second							
			TOTAL K	G / ELONG	ATION %				
COLOUR	SIZE IN MM	100%	200%	300%	400%	% 500% 1 29,2 0 66,3 2 97,0			
Peach	1040 x 15 x 3,0	7,3	10,0	14,4	19,1	29,2			
Orange	1040 x 22 x 4,0	13,0	20,7	27,1	35,0	66,3			
Lime Green	1040 x 30 x 4,4	22,5	33,0	44,0	56,2	97,0			
Blueberry	1040 x 40 x 4,4	25,6	40,0	53,0	76,0	132,0			
Plum	1040 x 53 x 4,4	38,0	56,0	79,0	105,0	167,0			

With the MoVeS F!T Superloop you can easily stretch and strengthen different muscle groups. It can be used for rehabilitation, fitness and strength training. The lay-flat-length from the MoVeS F!T Superloop is 104cm long (circumference is 208cm) and has an incredible elasticity.

The MoVeS FIT Superloops are multi-layered, this results in a more durable product. They are available in 5 different resistances: Peach, Orange, Lime Green, Blueberry and Plum.

			1. 1. 1.
COLOUR	RESISTANCE	MOVES F!T SUPERLOOP	UNIT
Peach	Light	01-110402	1 pc
Orange	Medium	01-110403	1 pc
Lime Green	Heavy	01-110404	1 pc
Blueberry	Extra Heavy	01-110405	1 pc
Plum	Special Heavy	01-110406	1 pc
			1991 2 4 4 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1

18

Our high quality and cost effective loops are also available in F!T colours. MoVeS F!T Loops are easy to use and very versatile.

Our loops are multi-layered, this results in a more durable product. Our thicker loops won't roll-up during exercising, avoiding the uncomfortable feeling from thinner alternatives.

Our loop set combines our most popular colours in a compact mesh bag. This affordable set is very useful to bring

Multi-Layered

pcs



MoVeS F!T Loop MD (E MoVeS F!T

01-RESISTIVE








MoVeS





COLOUR	RESISTANCE	MULTI-GRIP EXERCISER 30 Yard Roll (27,4m)
Tan	XX-Light	01-703001
Yellow	X-Light	01-703002
Red	Light	01-703003
Green	Medium	01-703004
Blue	Heavy	01-703005
Black	X-Heavy	01-703006
Silver	XX- Heavy	01-703007
Gold	XXX-Heavy	01-703008





COLOUR	RESISTANCE	MOVES	UNIT
COLOOK	RESISTANCE	ANKLECISER	UNIT
Yellow	Light	01-500102	5 pcs
Red	Medium	01-500103	5 pcs
Green	Heavy	01-500104	5 pcs
Blue	Extra Heavy	01-500105	5 pcs
Set of 4 resistances		01-500111	4 pcs

The MoVeS Bar is the ideal product to treat and resolve golfer's and tennis elbow. This effective product is developed to focus on hand, arm and shoulder problems, this product allows flexion and extension exercises, adduction and abduction exercises and even soft tissue mobilisation, far more than the twisting and bending you would imagine.

The MoVeS Bar is available in 6 colours / resistances: extra light, light, medium, heavy, extra heavy and special heavy.

> Ideal for Golfer's &

Tennis Elbow treatment

Size: 31cm x 4,5cm Packaging: polybag with full-colour header Material: dried natural rubber

1 size Ø; 6 Colours / Resistances

COLOUR	RESISTANCE	KG-FORCE*	MOVES BAR
Tan	Extra Light	2,5 - 3,5	01-600101
Yellow	Light	3,5 - 4,5	01-600102
Red	Medium	4,7 - 5,7	01-600103
Green	Heavy	6,7 - 7,7	01-600104
Blue	Extra Heavy	9,8 - 10,8	01-600105
Black	Special Heavy	12,0 - 14,0	01-600106
a angled other handles to a U share.			

*force applied when bending to a U-shape

The MoVeS Cuff-Ring Tube is a versatile exercise tool to tone and strengthen all the muscles of the lower body. With velcro closure that fits most sizes.



For intensive training of isolated muscles. Efficient for all muscles of upper- and lower body. The centre of the MoVeS 8-Ring Tube offers lots of exercises.

Ideal for abdomen, bottom, hip and leg exercises. With the MoVeS O-Ring Tube you get an intensive training of each muscular group.

COLOUR	RESISTANCE	O-RING	FIGURE 8-RING	CUFF-RING
	RESISTANCE	TUBE LOOP	TUBE LOOP	TUBE LOOP*
Yellow	Light	01-200202	01-200302	01-200402
Red	Medium	01-200203	01-200303	01-200403
Green	Heavy	01-200204	01-200304	01-200404
Blue	Extra Heavy	01-200205	01-200305	01-200405
Black	Special Heavy	01-200206	01-200306	01-200406





08-EVALUATION

01-RESISTIVE

02

-HAND

06-STRENGTH

07-THERMAI

MoVeS Ankleciser <€

Gold



14,7


The MoVeS Ankleciser is a great ankle exerciser, complying to all the rigorous quality standards MVS In Motion demands of its MoVeS products. It's as easy in use as it's effective. The feet slip comfortably into the latex free elastic band which consists of two stitched loops. Available in 4 different resistances, this progressive exercise system enables the user to choose the desired resistance level to perform the exercises.

5,9

Available in a set of 4 resistances, or in a 5-pack of the same resistance. Packed in a polybag with full-colour header and exercise chart.



01-RESISTIVE



02

22

01-RESISTIVE

04-WELLNESS

02-HAND



This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and notch handles.

The unique notch handles allow the user to easily adjust the length of the rope through the handle while the revolutionary stopper system is a great aid for consistent and progressive exercising.

The VALUE PACK is packed in a polybag, with notch handles but without stoppers on the nylon cord.

This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and wooden handles.

This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and detachable hand support.

This hand support allows patients without gripping power to use the shoulder pulley.

	Unit
h Handles Stoppers	1 pc
E PACK (notch handles, polybag, no stoppers)	1 pc
den Handles Stoppers	1 pc
Support Stoppers	1 pc



MoVeS Accessories \in





The spring loaded buttons of the MoVeS Flex-Ion isolate each digit, so the strong fingers cannot compensate for the weaker ones. This results in a better, faster and easier to control rehab process. MoVeS Flex-Ion is available in 6 colour coded units with gradual resistance. This gradual resistance helps the therapist measure progress. The ergonomic base adds comfort and stability. Comes with a full colour and multilingual instruction manual.



set

1 pc

MoVeS Flex-Ion Set | 6 Resistances | 1 Display 02-000111 MoVeS Flex-Ion Empty Display | Metal 02-000121









Reference	Description	Unit
01-300001	MoVeS Wall Mount (Plug & Play)	1 pc
01-300101	MoVeS Safety Handles	pair
01-300102	MoVeS Soft Handles	pair
01-300103	MoVeS Door Anchor	1 pc
01-300104	MoVeS Universal Assist	1 pc
01-300106	MoVeS Band & Tubing scissors	1 pc







01-RESISTIVE

04-WELLNESS

MoVeS Flex-Ion $\mathbb{MD} \subset \in$



08-EVALUATION

04-WELLNESS

05-FITNESS

06-STRENGTH

07-THERMAL



exercise leaflet.

COLOUR

Yellow

Red

Green

Blue Black Digi-Flex MD C€



The Handmaster Plus[™] strengthens the 9 muscles that close the hand as well as the 9 muscles that open the hand. The result is maximum strength, and balance to the hand, wrist and elbow.



Reference Description

02-050102	Handmaster Plus Soft Purple Ball - Yellow Cord
02-050103	Handmaster Plus Medium Red Ball - Black Cord
02-050104	Handmaster Plus Firm Orange Ball - Grey Cord
02-050111	Handmaster Plus Set of 3 (one of each level)

Constant Force X-Tend CE

Unlike elastic-based exercisers that increase in resistance the further the unit is extended, Constant Force X-Tend uses retractable-based resistance, which keeps a constant level of resistance from start to finish.

Ergonomic in design, the Constant Force X-Tend works the Intrinsic and Extrinsic Extensors of the hand and can isolate MP/PIP/DIP joint motion.

It also allows each finger to exercise independently with smooth, consistent resistance, thus improving range of motion, endurance, strength and dexterity more effectively.

- intrinsic and extrinsic extensor muscles throughout their full range of motion.
- well as the hand, wrist and forearm.

COLOUR	RESISTANCE (per digit)	REFERENCE	UNIT
Light	0,11 kg	02-100403	1 pc
Medium	0,18 kg	02-100404	1 pc
Heavy	0,27 kg	02-100405	1 pc



Digi-Flex Set | 5 Resistances | 1 Display

Digi-Flex Empty Display | Metal

Digi-Extend MD (€

REFERENCE UNIT

02-010103 1 pc 02-010104 1 pc

02-010105 1 pc

02-010121 1 pc

1 pc

1 pc

set

02-010102

02-010106

02-010111

The Digi-Extend[®] is a unique therapy device that can isolate the extension of all five fingers and each of the individual finger joints. The product strengthens the intrinsic and extrinsic muscle groups to maximise extensor tendon gliding, joint motion, coordination and muscle function in the hand and forearm.

Digi-Flex[®] develops isolated finger strength, flexibility and coordination.

Suitable for rehabilitation after hand injury. The Digi-Flex[®] comes with an

RESISTANCE (1 digit | all digits)

0,7kg | 2,3kg

1,4kg | 4,5kg

2,3kg | 7,3kg 3,2kg | 10,4kg

4,1kg | 14,1kg

The Digi-Extend® focuses on the extensor mechanism of one or all fingers as it targets the area of the hand that needs rehabilitation. Patients can easily change the resistance levels and vary their hand positions for specific exercises to help maximise the exercise goals and prevent injury. The Digi-Extend[®] can also be used for finger adduction, abduction, blocking and flexion exercises. It is compact in design, hightech, easy and comfortable to use as well as effective in helping patients reach their therapy goals. Four colour coded band resistance levels allow each muscle or joint to work at a comfortable level and build as increased strength is achieved.

Reference	Description	Unit
02-010201	Digi-Extend Bands (Tan to Green) included	1 pc
02-010211	Digi-Extend Set 4 units 1 Display	set
02-010202	Digi-Extend Band Set 25 pcs 5x Tan, 5x Yellow, 5x Red, 5x Green, 5x Blue	set
02-010221	Digi-Extend Empty Display Metal	1 pc





por m rex

VIA[™] is an entry level hand exerciser designed with young athletes in mind. It has a round, soft look, wide finger pistons for smoother function & improved finger performance.

The VIA[™] is available in two tensions.

and the second			
RESISTANCE	RESISTANCE (per digit)	VIA	UNIT
Light	1,5 kg	02-060102	1 pc
Medium	2,5 kg	02-060103	1 pc

Prohands Series $\zeta \in$



PRO[™] provides serious athletes with the most advanced hand exerciser in the world. The PRO[™] features a larger profile, superior steel springs that resist wear, an integrated ergonomic palm bar for comfort.

The PRO[™] is available in three tensions.

and the second se			
RESISTANCE	RESISTANCE (per digit)	PRO	UNIT
Medium	3,5kg	02-060303	1 pc
Heavy	4,5kg	02-060304	1 pc
Extra Heavy	5,5kg	02-060305	1 pc

_	
\leq	
in	
in	

30

01-RESISTIVE

02-HAND

03-UPPER & LOWER BODY

04-WELLNESS

08-EVALUATION

26

Handmaster Plus



CE



Unit 1 pc 1 pc 1 pc Set

Exercise each finger or thumb independently with smooth, consistent resistance to the

Improve range of motion, strength, dexterity, endurance, circulation for each finger as

Ergonomic design fits comfortably in hand, and includes ability to isolate pip & dip joints.





03-UPPER & LOWER

01-RESISTIVE

02-HAND

Varigrip Therapy $\in \in$

Resistance adjustable per finger

po m ex

The VariGrip Therapy has the ability to adjust tension in each individual finger, the resistance from each single finger knob can easily be adjusted. Resistance can be increased incrementally on the patient's need, which minimises overcompensation and strain. The VariGrip allows the patient to progress steadily between resistance levels.

5 different resistance levels are available.

Comparison Chart: VariGrip vs Digi-Flex [®] / Flex-Ion (kg per digit)										
VariGrip	X-Li	ght	Lig	ght	Mec	dium	He	avy	X-He	eavy
Digi-Flex [®] / Flex-Ion	0,7	1,4	1,4	2,3	2,3	3,2	3,2	4,1	4,1	n/a

Each Varigrip's highest resistance ends where the next VariGrip resistance starts. Any other hand exerciser increases from Light (1,4kg) to Medium (2,3kg). That's almost a 70% increase gap, which can cause an overload on the muscles, resulting in injuries. The seamless progress from one VariGrip to the other, supported by the padded finger and palm area, guarantees a safe recovery and better feel for those with sensitivity issues.

COLOUR	RESISTANCE (range per digit)	LEVEL	REFERENCE	UNIT
X-Light	0,7kg - 1,4kg	1	02-100102	1 pc
Light	1,4kg - 2,3kg	2	02-100103	1 pc
Medium	2,3kg - 3,2kg	3	02-100104	1 pc
Heavy	3,2kg - 4,1kg	4	02-100105	1 pc
X-Heavy	4,1kg - 5,4kg	5	02-100106	1 pc
VariGrip Set	5 Resistances 1 Disp	olay	02-100107	set
VariGrip Emp	02-100108	1 pc		



Varigrip Sport <€





COLOUR /	RESISTANCE	POWER-WEB COMBO	The F
Tan/Ultra Light	Red/Medium	02-020101	flexib
Yellow/Light	Green/Heavy	02-020102	to th
Red/Medium	Blue/Super Heavy	02-020103	const
Green/Heavy	Black/Ultra Heavy	02-020104	_





thumb.





effectively for individual fingers, as well as the entire hand, wrist and forearm.

VariGrip Sport has the patented ability to adjust the tension level in each individual finger, allowing the user to create a custom programme to build and improve endurance, strength and dexterity more

a callus builder, ideal for climbers, gymnasts, gamers, musi-

02-100304

and the second se			
Reference	Description	Resistance (kg)	Unit
02-100303	VariGrip Sport	2,3 - 4,1	1 pc
02-100304	VariGrip Sport Pro	4,1 - 5,4	1 pc

Both flexion and

Description

02-100204 VariGrip Uno

Reference

extension

Varigrip Uno <€

The only thumb/individual finger exerciser with adjustable resistance for both flexion and extension. An ergonomic design with stabilisation points offers extra support.

Gradually adjust resistance and better track progression based on the patients needs. Benefits are improved strength, dexterity, endurance and

range of motion.	Light Spring	0,7kg - 2,3kg		
	Heavy Spring	2,3kg - 3,6kg		
rip Uno includes two easily interchangeable springs.				

The VariGrip Uno The red spring allows a light resistance for early therapy or light training. Once a higher resistance is required the red spring can be swapped by the green spring, allow

U


1

	equired the red spring can
wing	a higher resistance.
nit	
рс	

Adjustable resistance allows a smooth transition between resistance levels



02-100303



02-HAND

07-THERMAL



28

Power-Web Combo MD C€



The #1 hand exercise web in the world!

Each Power-Web comes with a 39 exercise poster





Power-Web[®] Combo features a unique advantage two resistance levels in one, offering greater bility and value. It is identical in size and weight ne Power-Web[®] Senior with the same high quality struction and materials.

Power-Web Flex-Grip MD ⊂ €

The creators of the best-selling Power-Web[®], are proud to present the original Power-Web® Flex-Grip®. The Power-Web[®] Flex-Grip[®] is a unique exerciser that fits in your palm and provides the ability to work dexterity, extension and strengthening in your hands, fingers, and thumbs. Available in 3 different resistances, the Power-Web® Flex-Grip[®] has finger and thumb holes and has two sides: one side featuring a ball to squeeze and the other featuring buttons and a thumb pad to exercise fingers separately and increase dexterity and strength in your fingers and

Latex-free and easy to use!

Power-Web Senior $\mathbb{MD} \subset \in$

The Power-Web[®] Senior is 36cm in diameter and weighs 450g. Its unique patented design allows a variety of exercises to be performed comfortably and effectively. Constructed of high quality rubber with special agents added for durability and strength, the Power-Web® will last for years.

Now also available as latex-free!



MoVeS Squeeze Balls & Eggs MD CE



The MoVeS Squeeze Ball is developed to provide a progressive sequence of 5 colours. It strengthens hands, fingers, forearms and will work the muscles separately or in tandem. The MoVeS Squeeze Ball can be used for rehab as well as for sports training or stress relief.

The MoVeS Squeeze Egg is ergonomically shaped to provide highly effective exercise for the rehabilitation and strengthening of weakened or injured fingers, hands and wrists. Because it is small and easy to use, you can exercise anywhere, anytime.

Available colours/resistances: yellow, red, green, blue and black.

- The MoVeS Squeeze Balls & Eggs can help to:
- Rehabilitate post-operative/post-traumatic hands and arms.
- Restore function in hands, arms and shoulders following stroke.
- Reduce pain and stiffness in arthritic fingers, hands and wrists.
- Relax tense, sore muscles and reduce stress.
- Strengthen aging hands and arms.

cb m ex

COLOUR	RESISTANCE	LEVEL	SQUEEZE BALL
Yellow	Extra Soft	1	02-030102
Red	Soft	2	02-030103
Green	Medium	3	02-030104
Blue	Firm	4	02-030105
Black	Extra Firm	5	02-030106



COLOUR

Tan

Yellow

Red

Green

Blue

Black

RESISTANCE

Ultra Light

Light

Medium

Heavy

Super Heavy

Ultra Heavy

The Power-Web[®] Junior is 18cm in diameter and

weighs 225g. Just half the size of the Power-Web®

Senior, it has all the same features as the larger product and delivers many of the same benefits when used in a hand and wrist exercise programme.

cb m ex

01-RESISTIVE

02-HAND

PER & LOWER BODY

04-WELLNESS

POWER-WEB

SENIOR

02-020201

02-020202

02-020203

02-020204

02-020205

02-020206

POWER-WEB

LATEX-FREE SENIOR

02-022201

02-022202

02-022203

02-022204

02-022205

02-022206

Power-Web Junior MD C€



Adjustable resistance

Reference	Description	Unit
2-070201	MoVeS Wrist Exerciser	1 pc





MoVeS Wrist Exerciser <€


02-030205 02-030206



The MoVeS Wrist Exerciser is an adjustable wrist and forearm trainer which allows the user to exercise with progressive resistance. Its ergonomically designed forearm support keeps the wrist in the right place. Together with the soft handgrip, adjustable to 4 positions, it ensures maximum comfort while exercising your wrist and forearm.

02-HAND

01-RESISTIVE





22-88lb 10-40kg



The MoVeS Comfort Putty does everything a great putty should do: it strengthens weakened grasp through exercise and is used for rehabilitative exercise or stress relief. It is odourless, does not stain and will not leave any residue on the patient's hands. MoVeS Comfort Putty is non-toxic, greaseless, hypo-allergenic and gluten and latexfree. MoVeS Comfort Putty is available in 6 colours/ resistances and comes pre-packed.








COLOUR	RESISTANCE	IFVFI			THERAFL	EX PUTTY			
COLOUR	R RESISTANCE	RESISTANCE	LEVEL	57g	85g	113g	170g	454g	2,3kg
Tan	Extra Soft	0	02-040001	02-040101	02-040201	02-040301	02-040401	02-040501	
Yellow	Soft	1	02-040002	02-040102	02-040202	02-040302	02-040402	02-040502	
Red	Medium	2	02-040003	02-040103	02-040203	02-040303	02-040403	02-040503	
Green	Firm	3	02-040004	02-040104	02-040204	02-040304	02-040404	02-040504	
Blue	Extra Firm	4	02-040005	02-040105	02-040205	02-040305	02-040405	02-040505	

Y	e e andre		A r	djustable esistance
nce	Description		Unit	
102	Mambo Max Adjusta	able Hand Grip - Light	1 pc	





Foam Hand Grip

The Mambo Max Foam Hand Grip assures a firm training grip. The foam handles provide a firm grip for the entire workout session. Take it with you to the office or just use it at home. Wherever you want to use it, the Mambo Max Foam Hand Grip assures a tight grasp.



Reference Description

04-WELLNESS

01-RESISTIVE

02-HAND

- Why Anti-Microbial Agents? -

Unit

1 pair

Anti-microbial agents kill bacterium without additives and resist growth of odour-causing microbes. MoVeS Comfort Putty contains natural anti-microbial without addition of toxic chemicals. This additive ensures the MoVeS Comfort Putty material retains the clean, non-greasy feel from day one. Cheaper alternatives on the market do not contain these agents. As a small amount of bacteria and other micro-organisms will always be present on the patient's hands and because putty should be stored in room temperature conditions, anti-microbial agents are a must to avoid the putty becoming the ideal breeding ground for bacteria and fungi. Anti-microbial agents extend the product's lifespan and improves its cleanliness.







T PUTTY
454g
02-090401
02-090402
02-090403
02-090404
02-090405
02-090406







Theraflex Putty MD C€

Theraflex[™] Putty strengthens weakened grasp through exercise and is used for rehabilitative exercise or stress relief. It is odourless, does not stain and will not leave any residue on the patient's hands. Theraflex[™] Putty is nontoxic, greaseless, hypo-allergenic and latex-free.

Theraflex[™] Putty is available in 5 colours/resistances.



Fit Stretch **MD** ⊂ €





III Cb

The Mambo Max Fit Stretch is a fantastic, simple and portable tool to stretch your lower leg muscles. It is perfect for warm-ups and cooldowns. The use of a Mambo Max Fit Stretch helps increase your ankle, shin and calf flexibility. It is made of durable, moulded plastic and includes a non-skid foot pad.

Reference	Description	Unit
03-050101	Mambo Max Fit Stretch Single Leg	1 pc



Because the Saehan Cycle - Clinical

is a product aimed at professionals,

we offer various accessories to

increase its flexibility and user's safety. Be aware these do only apply to the clinical cycle, not to

the home cycle.

Saehan Cycle (**SAEHAN**

The Saehan Cycles are pedal exercisers aimed at increasing strength and flexibility.

Both the Home and Clinical Saehan Cycles include attachments to fit the pedal exerciser to a chair or even a wheelchair. The pedal exerciser's length can be perfectly adjusted to the user when his feet are fitted in the pedals with foot straps. A convenient turning knob makes it possible to adjust the resistance (0 - 13 kg) which increases or decreases the effort, resulting in a progressive exercise programme.

The Clinical version from the Saehan Cycle has an even higher level of adjustability. The distance between pedal and crank can be adjusted and the resistance knob indicates the level selected.

Reference	Description	Unit
03-010311	Clinical Cycle Accessory Long front leg	1 pc
03-010312	Clinical Cycle Accessory Extension arm (right = left)	1 pc
03-010313	Clinical Cycle Accessory Hand strap - right	1 pc
03-010314	Clinical Cycle Accessory Hand strap - left	1 pc
03-010315	Clinical Cycle Accessory Long (straight) leg brace set	1 pc
03-010316	Clinical Cycle Accessory Bent leg brace	1 pc

Reference	Description	Unit
03-010301	Saehan Cycle - Home	1 pc
03-010302	Saehan Cycle - Clinical	1 pc





Wheeleo MD (E SWheeleo Yesterday, you had the choice between a classic cane which is impractical because it is difficult to coordinate, does not help when it is in the air. No one likes having to use it because it highlights their difficulties... Or a rollator which requires two hands to use it, which stigmatizes as much as the cane, which is cumbersome and not easy to handle, let alone indoors (try walking through a door with a rollator!). Therefore, the walker is parked at the entrance of the house or apartment and leaves its user without help. Now there is the Wheeleo.

Reduced mobility presents 3 challenges:

- Keep the person moving; •
- Reduce the risk of falling;
- Loss of autonomy.

Wheeleo, the rolling cane, is designed to offer a solution to these 3 challenges and offers additional benefits:

- Ease of handling, no special training required;
- Walking goes a lot faster than with a quadripod cane;
- Safer walking than with a quadripod cane, as there is constant ground contact;
- One-handed use while still maintaining constant ground contact, the other hand can be used to move an object, to open doors, ...
- Minimal space required and perfect for indoor use where a rollator is too big.
- Reinforced stability: the wheels offer a security system. A lateral push will result in a lateral movement, not a tilting/fall.

user.



Rolling Cane

Wheeleo

Adjustable height

> \$ 77,5 - 101,5cm ↔ 31 x 27cm 120kg

possibly fall.

lifting the Wheeleo.

02

03-UPPER & LOWER BODY

04

0





03-UPPER & LOWER BODY



EVALUAT -00



Wheeleo is a unique quadripod rolling cane used to improve mobility. It is designed with the aim of offering help without stigmatising the disability of its

It will make walking trips smoother and safer which will allow the user to maintain or improve autonomy.

The Wheeleo is used with one hand. Thanks to its 4 wheels, you do not have to lift it, which improves stability and reduces the risk of falling.

The Wheeleo is super easy to use. Ultra-intuitive, its user will always position it in the right place to guarantee its stability and balance.

An important feature of the Wheeleo is that it **does not require brakes** to prevent its user from falling. Thanks to its rolling design, the Wheeleo will always be placed in the right position for its user.

When losing balance, the user will automatically push the Wheeleo in the right position that is

comfortable for him, without

Adding brakes would be the same as replacing the Wheeleo with a regular cane because the brakes would block the wheels. As a consequence the Wheeleo would tilt and its user would



MoVeS Oxycycle 1 - Active Magnetic MD CE



M_oV_eS

The OxyCycle 2 comes with

half-foot pedals and


handgrip pedals.

MoVeS Oxycycle 3 - Active & Passive <€





Description	Unit
MoVeS Strap Pedals	pair
MoVeS Half-Foot Pedals	pair
MoVeS Full-Foot Pedals	pair
MoVeS Handgrip Pedals	pair
	MoVeS Strap Pedals MoVeS Half-Foot Pedals MoVeS Full-Foot Pedals



For PASSIVE training only

Adjustable

ASSISTANCE

The OxyCycle 1 Pedal Exerciser is used for upper and lower body exercises. The adjustable resistance makes it a great product to progressively strengthen your arms and legs. The OxyCycle 1 Pedal Exerciser has a light-weight compact design, adjustable resistance and a multi-function LCD display readout of scan, time, distance, calories and RPM. A non-skid mat and interchangeable handgrip pedals are included.

The MoVeS OxyCycle 2 Pedal Exerciser is a great tool for personalised, low-stress workouts at home. The adjustable speed level makes it a great product to progressively

strengthen your arms and legs. This pedal exerciser

provides low impact exercise to increase strength and

flexibility, is versatile for both upper and lower body and it improves the blood circulation. It offers up to 60 Watt of

motorized assistance. It has enough power to maintain a constant rotational speed, whereas cheaper models only carry a 30 Watt motor. A safety stop will avoid injury in case of cramp or muscle spasm. The MoVeS OxyCycle 2 has a light-weight compact design, adjustable speed and a multi-

function LCD display readout of scan, time, count, RPMlevel and estimated calories burned. A non-skid mat and

interchangeable handgrip pedals are included.







For home use only

For home use only





• 50x46x38 cm Comes with UK or EU plug!

- 8,5 kg
- 50 Hz
 - 25 65 rpm

Reference	Description	Unit
03-010102	MoVeS OxyCycle 2 Pedal Exerciser EU-plug	1 pc
03-010102 UK	MoVeS OxyCycle 2 Pedal Exerciser UK-plug	1 pc

HINESS

02

03-UPPER & LOWER BODY

04

36





For ACTIVE and PASSIVE training

The MoVeS OxyCycle 3 Pedal Exerciser is a great tool for personalised, low-stress workouts at home. The adjustable speed level and resistance makes it a great product to progressively strengthen your arms and legs. This pedal exerciser provides low impact exercise to increase strength and flexibility, is versatile for both upper and lower body. It has a reliable **80 Watt** motor to offer enough power to maintain a constant rotational speed, whereas cheaper models only carry a 30 Watt motor. A safety stop will avoid injury in case of cramp or muscle spasm. The MoVeS OxyCycle 3 has a light-weight compact design, adjustable speed & resistance, a 15-minute timer switch and a multifunction LCD display readout of scan, time, count, total count and calories burned. A non-skid mat and interchangeable

MoVeS Oxycycle Pedals $\in \in$







Half-Foot Pedal

Curamotion *C*€

W Cb m EXER 1 96 x 45 x 45 (cm) / 11,5 kg

Wall mount included

146 x 45 x 54 (cm) / 13,5 kg Wall mount included Adjustable height

EXER 2

Improvement of muscle strength and joint range of motion. The CuraMotion Exer Series has an adjustable resistance and a multi-functional LCD-display readout scan which indicates rotation count, rotations per minute, time exercised and calories consumed.

All are easy adjustable with a quick release knob.





EXER 3 37 x 38 x 70 (cm) / 9 kg With foot stands

Reference	Description	Unit
03-030101	CuraMotion Exer 1	1 pc
03-030102	CuraMotion Exer 2	1 pc
03-030103	CuraMotion Exer 3	1 pc
03-030104	CuraMotion Exer 4	1 pc



EXER 4 96 x 45 x 45 (cm) / 13 kg Wall mount included



MoVeS Pedal Exercisers are ideal for restoring impaired muscle strength, improving blood circulation and give a better coordination to the legs or arms. It can be used on a table top as an arm exerciser or on the floor to give the legs a workout. A convenient resistance control knob allows a broad range of pedal resistance for a progressive exercise programme.

- ٠ and calories)
- •



MoVeS Standard Pedal Exerciser | Premium Chrome

- Pre-assembled •
- Pedal straps adjustable to two positions ٠
- Quality chrome finish ٠
- In colour box or on specific demand in bulk with separate colour box ٠



MoVeS Standard Pedal Exerciser K/D

- Knock-down for compact shipping and storage (only 6 screws to assemble or disassemble) •
- Pedal straps adjustable to two positions
- Budget solution ٠



0

80

02 HAND

_NESS

MoVeS Pedal Exerciser MD C€




01-RE

03-010203



Unit	
1	рс
1	рс
1	рс



MoVeS Wooden Incline Board $MD \in E$

The incline board increases extension in the lower and upper back while stretching the calves. It can be used:

• to improve the general posture,

///

- to stretch the calf muscles of athletes, football players and sportspersons in general, •
- to stretch the upper-leg muscles (Ischias),
- to increase extension in the lower and upper back.

MoVeS Wooden Incline Board - Ladder



The MoVeS Wooden Incline Board - Ladder is made out of natural birch wood and is very resistant. It is adjustable in 4

different angles (10°, 15°, 20° and 25°). After using the MoVeS Wooden Incline board, it can easily be folded to achieve its compact form when not in use

MoVeS Wooden Incline Board - Fix



The MoVeS Wooden Incline Board - Fix is perfect for who's looking for an easy and budget-friendly Incline Board. It's made of polymerised wood for a strong feel and high user weight. The inclination (17,5°) is perfect for physiotherapy exercises.

Reference Description	Unit
03-040001 MoVeS Wooden Incline Board - Ladder	1 pc
03-040003 MoVeS Wooden Incline Board - Fix	1 pc

Plastic Incline Board MD C€

Unit

1 pc

Apart from stretching the calf and upper-leg muscles, the incline board is a fantastic device to improve general posture and extension in the lower back.

Made of a composite plastic for a maximum load of 250kg. Adjusts to 4 levels and folds flat (8cm) to store easily. Angles are 5°, 15°, 25° and 35°. The last slot on the board makes an excellent hand-hold for stretching. Non-slip top measures 35x35 cm.

՝ 250kg	
Reference	Description
03-040102	Incline Board Standard





Adjustable

50cm x 29,5cm x 3,5cm

0 150kg



Kinesiology Tape MD C€



MoVeS Kinesiology Tape is made from quality cotton fibres and is used for joint and muscle support. This skin-friendly Kinesiology Tape is hypoallergenic, latexfree, breathable and water resistant. The adhesive side of the Kinesiology Tape is made of acrylic which prevents sweat and moisture from penetrating the tape. This tape with a wave-like structure can be stretched up to 170% lengthways.

Available in 4 colours.



Reference	Description		Unit
04-060101	MoVeS Kinesiology Tape 5cm x 5m Tan	NEW	1 pc
04-060102	MoVeS Kinesiology Tape 5cm x 5m Black	NEW	1 pc
04-060103	MoVeS Kinesiology Tape 5cm x 5m Pink	NEW	1 pc
04-060104	MoVeS Kinesiology Tape 5cm x 5m Light blue	NEW	1 pc

AQ[®] Massage Lotions SKIN IN BALANCE

NAQI[®] Massage Lotions

- Excellent Lubrication
- Economical usage
- Hypo allergen perfume Prolonged hydration
- Skin-friendly



02-HAND

07

THERMAI

80

40





NEW

YOUR HANDS DESERVE THE BEST POSSIBLE

Rinses off with water

Non greasy leaves skin soft and drv

smooth, ight skir feel

02-HAND

04-WELLNESS

NAQI[®] nv Stadsbeemd 1037, 3545 Halen, Belgium

www.nagi.com

Tel.+32(0)13.460.190





Mambo



Massage rollers are used to put pressure on the muscles' triggerpoints and to loosen the tension in the muscles. Trigger points are often the cause of different pains.

Spiky Foam Roller



Hollow Foam Roller

The Mambo Max Hollow Foam Roller can easily help people who are experiencing long term pains. Its three-dimensional pattern with different shapes replicates the massage from a therapist's hands. The larger surfaces allow a superficial massage while the smaller shapes are designed for a deeper massage.



Ergonomic Foam Roller



3-in-1 Foam roller

The Mambo Max 3-in-1 foam roller combines 3 different foam rollers in one package. Thanks to its smooth surface, the outer part is perfect for a superficial massage. The second part has grooves which can be used for a deeper massage and to exercise higher pressure on the trigger points. Finally, with the inner massage stick it's possible to exercise high pressure on sore muscles to improve blood circulation and the muscles' recovery time. This last item is available as a separate product too: The Mambo Max Massage Stick.

Reference	Description	Unit
04-050101	Mambo Max Hollow Foam Roller	1 pc
04-050102	Mambo Max Ergonomic Foam Roller	1 pc
04-050103	Mambo Max Spiky Foam Roller	1 pc
04-050301	Mambo Max 3-in-1 Foam Roller	1 set
04-050201	Mambo Max Massage Stick	1 pc

	- 1 yoga strap with buckle (180 x 3,8 cm)
	This set is packed in a full colour box. Items are also sold separately.
Mambo Nga Set View View	

Reference	Description	Unit
04-010201	Mambo Max Yoga Mat	1 pc
04-010202	Mambo Max Yoga Block	5-pack
04-010203	Mambo Max Yoga Strap	1 pc
04-010211	Mambo Max Yoga Set 1 Mat 2 Blocks 1 Strap	Set

Stepping Stones $\in \in$

The non-slip 16 cm Mambo Max Stepping Stone is an amazing tool to spice up your workout. They can be used by both children and adults for indoor and outdoor workouts to improve dynamic balance, body awareness and coordination. It can also serve as a foot massage tool while seated.

Unit



The Stepping Stones now come in bright and vivid colours: yellow, orange, red, purple, blue and green.

They are packed in cardboard boxes per 6 pieces and each piece has a different colour.

ີ່ 100kg		
Reference	Description	

6-pack 04-030141 Mambo Max Stepping Stone


02

III ex

80

42

٠





The Mambo Max Spiky Foam Roller is designed for deeper massage & acupressure therapy. The spikes make it possible to get deeper into the muscles to relieve tension and soreness and to improve mobility. It also has a bi-directional zone for comfort when rolling along the spine. The Spiky Foam Roller is mainly used for stretching, self-massage, body support and stability exercises.

→ 33cm Ø 14cm



The Mambo Max Ergonomic Foam Roller is a curved hollow foam roller, which enables more exercises. Except rolling back and forth, it's possible to make sideways movements. More muscles will be addressed through which the body will be toned even more. The soft foam in combination with the nipples is perfect for trigger point massages. The central groove will reduce pressure on the spine.



01-RESI

02-HAND

LOWER BODY

06-STRENGTH

TRE

08-EVAL

43

Sitting Cushions <€



Ø 33cm

\$ 7cm

Can be used as

Balance Trainer

200kg


Ø 32cm \$ 7cm

200kg


Mambo Max Sitting Cushions are air-filled cushions for active and dynamic sitting, recommended by therapists. These easy-to-carry air cushions activate intrinsic trunk muscles to support the spine and are suitable for everybody. The surface from these sitting cushions can easily be cleaned with soap and water.

Comfort Cushion 33

111 CD m p



Standard Cushion 37



Faster Blaster pump **Sitting Wedge**



one smooth | one senso side

Reference	Description	Unit
04-020301	Mambo Max Comfort Cushion 33 Blue + Needle Pump	1 pc
04-020302	Mambo Max Comfort Cushion 33 Red + Needle Pump	1 pc
04-020303	Mambo Max Comfort Cushion 33 Black + Needle Pump	1 pc
04-021301	Mambo Max Comfort Cushion 33 Blue in Polybag	1 pc
04-021302	Mambo Max Comfort Cushion 33 Red in Polybag	1 pc
04-021303	Mambo Max Comfort Cushion 33 Black in Polybag	1 pc
04-020101	Mambo Max Standard Cushion 33 + Needle Pump	1 pc
04-020102	Mambo Max SensAir Cushion + Needle Pump	1 pc
04-020103	Mambo Max Standard Cushion 37 + Faster Blaster Pump	1 pc
04-020201	Mambo Max Sitting Wedge + Faster Blaster Pump	1 pc



Standard Cushion 33

SensAir Cushion

one smooth | one senso side

one smooth | one senso side

Soft-Over-Ball 👩 🛒 💇 The Mambo Max Soft-Over-Ball is an inflatable ball used for resistance work in Pilates. It is a 100kg

///_____

Reference	Description
04-010105	Mambo Max Pilates Soft-Over-Ball - 17-19 cm 0
04-010113	Mambo Max Pilates Soft-Over-Ball - 17-19 cm I
04-010107	Mambo Max Pilates Soft-Over-Ball - 17-19 cm I
04-010115	Mambo Max Pilates Soft-Over-Ball - 21-23 cm 0
04-010123	Mambo Max Pilates Soft-Over-Ball - 21-23 cm I
04-010117	Mambo Max Pilates Soft-Over-Ball - 21-23 cm I
04-010125	Mambo Max Pilates Soft-Over-Ball - 25-27 cm 0
04-010103	Mambo Max Pilates Soft-Over-Ball - 25-27 cm I
04-010127	Mambo Max Pilates Soft-Over-Ball - 25-27 cm I

Foam Roller 🐠 🌌

The Mambo Max Foam Roller has been designed for Pilates training and spine exercises (back aches) and is made of EVA. Ideal for balance and stabilisation exercises. Foam rollers are used for stretch and strength exercises, as well as for massage and relaxation.



Reference	Description	Unit
04-010102	Mambo Max Pilates Foam Roller - EVA 90cm	1 pc
04-010104	Mambo Max Pilates Foam Roller - EVA 45cm	1 pc

Ring of mex



Foam covered fiberglass resistive exercise ring with foam padded grips. For exercises focusing on strength, balance, flexibility and posture. The Mambo Max Pilates Ring adds light resistance to pilates exercises and helps target specific muscle groups. Often used to help tone the arms and thighs, and increase the core strength challenge of an exercise.

Ø 38 cm

Reference	Description	Unit
04-010101	Mambo Max Pilates Ring	1 pc



02-HAND

04-WELLNESS

44


Pilates <€



wonderful tool for adding focus and energy to your workouts. In addition, when slightly deflated,



it can be used as a soft prop under the head, neck, feet, back... Now available in 3 different colours and sizes, giving a range of 9 different Soft-Over-Balls.







07-THE TREAT 08-EVALUATION

45

Massage Tools <€

The Knobble® II massage tool is made of durable polymer with a band of non-latex. The Knobble® II gives your hand more power when you need it. Use either the tip or the base for applying pressure on yourself or others. The Knobble® II doesn't replace your healing hands, it just helps them do more and last longer.





The Jacknobber® II massage tool can be used in many ways to apply deep pressure to myofascial trigger points and areas of muscle tension. The smaller knobs provide more intense pressure while the larger knobs provide a more diffuse pressure. Use one knob at a time to press deeply into taut bands with embedded trigger points; use two knobs at a time on either side of the spine.

The Index Knobber[®] II massage tool is a simple device designed to let you apply deep pressure to muscles within easy reach of your hand. It is the ideal instrument for clinic or home use allowing for firm, sustained pressure on targeted muscle. It's unique design allows it to be used in several hand positions and by the right or left hands comfortably and interchangeably.





The Palmassager[™] is one of the most efficient massage therapy aids available anywhere. Hold it firmly in the palm of the hand and use one, two or three of its round knobs to massage sore muscles. It is used for back, shoulder, arms and legs.

pain-			
press	Unit	Description	Reference
as rec	1 pc	Jacknobber II	04-031101
achie	1 pc	Index Knobber II	04-031102
effect	1 pc	Palmassager	04-031103
conve	1 pc	Backnobber II	04-031104
_	1 pc	Knobble II	04-031105

The Backnobber® II massage tool is allows the user to apply as much -releasing sure as long

eauired to

ieve remarkable

cts. A locking mechanism allows it to come apart for enient travel and storage.

Massage Ball, Roll & Ring \in

The Mambo Max Massage Balls, Roll & Ring are used for relaxing muscles and reducing tension by relieving trigger points. The Mambo Max Massage Ball helps to strengthen hands and feet and to improve coordination and flexibility. It can be used for massage as well as to improve dexterity skills.



Reference	Description	Unit
04-030103	Mambo Max Massage Ball 6 cm	1 pc
04-030105	Mambo Max Massage Ball 7 cm	1 pc
04-030107	Mambo Max Massage Ball 8 cm	1 pc
04-030109	Mambo Max Massage Ball 9 cm	1 pc
04-030111	Mambo Max Massage Ball 10 cm	1 pc
04-030121	Mambo Max Massage Roll 16 cm	1 pc
04-030131	Mambo Max Massage Ring	1 pc



Mambo





Gym Mat



With rings

and carry strap.

The Mambo Max Gym Mats provide a soft touch and comfortable exercise environment. They are perfect for gym & fitness. Non-slip, well-cushioned surface with a unique closed cell structure. The mats are equipped with plastic rings 40cm apart so they can be suspended while not in use. Each mat comes with a handy storage/carry strap.

The Comfort Gym Mat comes in 3 popular colours: black, blue and magenta. Our well-known 1,5cm-thick exercise mats are now called 'Xtra Comfort Gym Mat' to highlight the extra comfort they offer compared to thinner mats. They continue to be available in widths of 60cm and 100cm in blue, red and green.



Nelerence	Description
05-030001	Mambo Max Gym Mat Wall Hanger for 10 mats
05-030303	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) I
05-030304	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) B
05-030305	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) I
05-030101	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (
05-030102	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (
05-030103	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (
05-030201	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5
05-030202	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5
05-030203	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5

///

01-RESISTIVE

02-

-HAND

///

46

07-THERMAL





Xtra Comfort Gym Mat 180 x 100 x 1,5 (cm) Xtra Comfort Xtra Comfor 1,5cm 1,5cm



Step aerobics is a form of aerobics that is practised worldwide by millions of

people. The reason for its popularity is simple, it offers a great cardio workout

While the aim of aerobics is to increase the body's consumption of oxygen,

without the need for extensive space or equipment.



Mambo



Fabric Elastic Band <€

The Mambo Max Fabric Elastic Band consists of 10 consecutive loops. This great toning and work-out tool is suitable for early training up to professional athletes. Three different resistances are available, while the 10 consecutive loops allow you to spice things up. Physical training can be intensified by using a lower amount of consecutive loops resulting in a higher resistance.



Mambo
Fabric
Elastic Band
77
10

total: 120cm loop: 11,5cm

COLOUR	RESISTANCE	REFERENCE	UNIT
Light	7 kg	05-120107	1 pc
Medium	10 kg	05-120110	1 pc
Heavy	15 kg	05-120115	1 pc

versatile product enables the user to perform exercises on both sides of the product, allowing it to be used for a full body work-out as well as for balance exercises. Basically it's Mambo Max' economical answer to the BOSU° system, but with extra tubing and handles!





step aerobics builds on that goal by offering a more intense workout. The Mambo Max Aerobic Step provides this extra boost by incorporating the action of stepping onto the platform to intensify the workout. It is low impact, and doesn't stress the joints as much as running or jogging. 15cm 20cm 25cm 96cm x 33cm 0 100kg \$ 15-20-25cm Reference Description Unit 05-070101 Mambo Max Aerobic Step 1 pc **Ball Accessories** CE /// Faster Blaster Pump Needle Pump Ball Stacker Electric Compressor Pump with velcro Hand Pump closure oot Pumr **Ball Stabiliser** Reference Description Unit 05-020101 Mambo Max Ball Stabiliser 1 pc The Ball Stabiliser facilitates balance. 05-020102 Mambo Max Ball Stacker 3 pcs Its notches enable combined AB 05-020201 Mambo Max Faster Blaster Pump | 2x120cc 1 pc Gym Ball / MoVeS Band exercises. 05-020202 Mambo Max Hand Pump | 2x2000cc 1 pc 05-020203 Mambo Max Needle Pump | 100cc 1 pc 05-020204 Mambo Max Foot Pump 1 pc

1 pc

01-RESISTIVE

111 Cb

BODY

07

05-020205 Electric Compressor Pump



The Mambo Max Fabric Elastic Band is durable and extremely versatile. It can be used for stretching, toning or even as an alternative for dumbbells as a weight training accessory. The progressive resistance offers you a broad range of exercises.





1 pc

08-EVALUATION

02-HAND

03-UPPER & LOWER BODY

06-STRENGTH

07-THERM

Balance Board <€



The Mambo Max Balance Board is an effective tool for balance training and fitness exercises. Stand on the rounded platform and start balancing. It is also a great tool to strengthen your ankles and thigh muscles. For a more severe balance training, add the included screw-on top to your Mambo Max Balance Board.

Ø 42cm \$ 5-8cm 0 100kg



Unit

1 pc

Description Mambo Max Balance Board

Wooden Balance Boards \in

The Mambo Max Wooden Balance Board is your perfect partner if you want to improve your balance on an ecologic and high-end product. The Wooden Balance board is available in 3 versions: Uni-Directional, Standard & Sport. Their different form and balance point make it possible to have 3 products with a very different character and level of difficulty.

Children and elderly could experience a balance board with pivot point as too challenging. The Uni-directional Balance Board offers more confidence. This enables the user to improve body stability without the fear of falling.

A Uni-directional Balance Board offers a more predictable training as it can only pivot in one direction. When turning the Balance board upside-down, the arches can be used as push-up grips. Once the user is ready for the next level of training, the Balance Board Rehab or Sport are ideal next steps in the process.

The Balance Board Standard is lower than the Balance Board Sport. In addition with its smoother pivot point, it reacts completely different than the Sport. It's a great middle step between the Balance Board Uni-Directional and Balance Board Sport.

The most challenging Balance Board is our Sport version, which, as the name specifies, makes it a real sport. Finding your balance and avoiding the edges touching the ground is hard and requires you to use all your muscles in your body. It's a great and fun way to exercise a fullbody workout.

Reference	Description	Unit
05-040001	Mambo Max Wooden Balance Board - Standard	1 pc
05-040002	Mambo Max Wooden Balance Board - Sport	1 pc
05-040003	Mambo Max Wooden Balance Board - Uni-directional	1 pc

Balance Board Freeman <€

1 pc

The Mambo Max Balance Board Freeman is an effective tool for balance training. Stand on the rounded platform and start balancing. This is a great way to strengthen your ankles and thigh muscles. The handgrips on the bottom are used for upper body balancing and workouts.

Ø 40c ≎ 8cn	ะm า		
Referenc	e Descriptio	on	
05-04020)2 Mambo M	Aax Balance Board Free	man

Uni: 40cm x 40cm / \$ 10cm 150kg



Standard: 40cm x 40cm / \$ 6cm Max inclination: 12,5° 0 150kg



Max inclination: 20° 150kg

Mambe.

Ecofriendly pvc. Heavy metal content in compliance with EN 71 Standards. Phthalates content meets EU-directive 2005/84/EC.

ANTI-BURST

350kg Yes, 350 kg!!!

Ball Size	Colour	Body Height	Mambo Max AB Gym Ball Box + Faster Blaster Pump	Mambo Max AB Gym Ball Polybag	Unit
45cm	Yellow	< 125 cm	05-010102	n/a	1 pc
55cm	Red	125-150 cm	05-010103	05-010203	1 pc
65cm	Green	150-175 cm	05-010104	05-010204	1 pc
75cm	Blue	175-185 cm	05-010105	05-010205	1 pc
85cm	Black	185-200 cm	05-010106	n/a	1 pc
95cm	Silver	> 200 cm	05-010107	n/a	1 pc

AB Peanut Ball <€



The Mambo Max AB (Anti-Burst) Peanut ball helps users who struggle with the instability of a regular gym ball. The peanut shape offers extra sideway stability, which gives more confidence to start working out with Gym Balls.





Peanut Ball Size Colour

Red

50 x 100cm

Mambo Max AB Peanu + Faster Blaster Pur 05-011103

02

Cb ex

05-040201

III ex

06-

RMAL

80



The Mambo Max AB (Anti-Burst) Gym Ball increases flexibility, improves strength, fitness and posture. Exercising with this ball also enhances your balance and coordination. All AB Gym Balls are Anti-Burst: if the ball gets punctured by a sharp object the ball won't burst, thus reducing risk of injury.

This extremely durable Mambo Max AB Gym Ball comes in a Mambo Max cardboard retail pack, including AB Gym Ball, Faster Blaster Pump, Exercise Chart, Plug and Plug Remover.



The Mambo Max AB Gym Ball is also available in the value pack: a polybag containing only the ball and plug.







The Mambo Max AB Peanut Ball offers a more predictable therapy or fitness training as it can only roll in one direction. This certainly is a huge benefit for children, elderly or people who suffered an injury. In addition, Peanut Balls are very popular amongst pregnant women and are very helpful devices during labour



Similar to our regular Mambo Max AB Gym Ball, the Peanut Ball is

Anti-Burst: if the ball gets punctured by a sharp object the ball won't burst, thus reducing risk of injury.

The Mambo Max AB Peanut Ball comes in a Mambo Max cardboard retail pack, including Faster Blaster Pump and Plug.

ıt Ball	Unit	
np	Onit	
	1 pc	

02-HAND

03-UPPER & LOWER BODY

NESS



02

03

07

52



	Unit
	1 pc
acite	1 pc
	pair

01-RESISTIVE



Coated Dumbbell < €

The Mambo Max Dumbbells are cast iron dumbbells, coated with vinyl or neoprene. The vinyl coating is easy washable or can be treated with a water based disinfectant.

The neoprene coating has an added resin to give them a scrubby feel and superior grip.

Available in 7 weights with matching colour code. The Dumbbell Rack holds 6 Dumbbells (3 pairs).



Description	Weight	Vinyl	Neoprene	Unit
Mambo Max Dumbbell	0,5kg	06-010102	06-010202	pair
Mambo Max Dumbbell	1,0kg	06-010103	06-010203	pair
Mambo Max Dumbbell	1,5kg	06-010104	06-010204	pair
Mambo Max Dumbbell	2,0kg	06-010105	06-010205	pair
Mambo Max Dumbbell	3,0kg	06-010107	06-010207	pair
Mambo Max Dumbbell	4,0kg	06-010108	06-010208	pair
Mambo Max Dumbbell	5,0kg	06-010109	06-010209	pair
Mambo Max Dumbbell rack (empty)	n/a	06-00	00101	1 pc

SoftMeds <€

111 Cbrex Mambo Max SoftMeds are colour coded hand weighted balls that are small and flexibleenough to grasp with one hand, yet big enough to hold with both hands. The soft, pliable vinyl cover makes it easy to hold and is a safe and effective way to add weight training to your workout.

Mambo Max SoftMeds are fully inflated and ready to use, straight out-of-the box. When needed, air pressure can be adjusted with a needle pump.



Ø 12cm

Reference	Weight	Description
06-040102	0,5kg	Mambo Max SoftMeds Tan
06-040103	1,0kg	Mambo Max SoftMeds Yellow
06-040104	1,5kg	Mambo Max SoftMeds Red
06-040105	2,0kg	Mambo Max SoftMeds Green

FAIRS Therapie Hamburg

Hamburg - Germany 22 - 23 October 2021

Therapie Hamburg is the platform in northern Germany for physiotherapists, occupational therapists and sports therapists

as well as speech therapists, masseurs and lifeguards. Innovations, products and industry trends are shown here.

Arab Health

therapie

ARAB HEALTH Dubai - UAE

24 - 27 January 2022

By Informa Markets Arab Health is the largest gathering of healthcare and trade professionals in the MENA region. With over 55,000 local and international visitors, 4,250 exhibitors and participants from over 150 countries, Arab Health is truly where the healthcare world comes to do business.

comed 5,100 exhibitors, nearly 130,000 visitors from 144 countries and had a satisfaction rating for closed deals of 96%.

MEDICA

GLOBAL FIBO FITNESS Cologne - Germany

FIBO Global Fitness is an expo for sports,

nutrition, and healthy lifestyle. It brings global top deci-

///

MEDICA

Düsseldorf - Germany

15-18 November 2021

MEDICA has been the largest

medical trade fair in the world

and the event on the agenda

of medical experts for over 40


7 - 10 April 2022

Mambo

years. In 2019, the fair wel-

sion makers of the industry together looking for the latest innovations and trends for their gyms, wellness facilities, doctors and physiotherapy practises, etc...

Swim Noodle $\in \in$



Ø 7cm ↔ 150cm

///



05-FITNESS


04-WELLNESS

080

54











Unit			
1	рс		





Wrist & Ankle Weight <€



Reference	Description	Weight	Unit
06-020102	Mambo Max Wrist & Ankle Weights	0,5kg	pair
06-020103	Mambo Max Wrist & Ankle Weights	1kg	pair
06-020104	Mambo Max Wrist & Ankle Weights	1,5kg	pair
06-020105	Mambo Max Wrist & Ankle Weights	2kg	pair
06-020107	Mambo Max Wrist & Ankle Weights	3kg	pair
06-020108	Mambo Max Wrist & Ankle Weights	4kg	pair
06-020109	Mambo Max Wrist & Ankle Weights	5kg	pair



To quickly build stamina and muscle, simply add Mambo Max Wrist and Ankle Weights to your workout. They are filled with iron sand and will effectively shape and tone your body and can be worn anywhere. Wear them around the office, out shopping or at the gym. Eventually, your body will adjust to running, walking and living with this extra weight.





Medicine Ball < €

A medicine ball (also known as an exercise ball, a med ball, or a fitness ball) is often used for strength training as it plays an important role in the field of sports medicine. Medicine ball training is one of the oldest forms of strength and conditioning training – the first reference to wrestlers training with sand filled bladders appears in Persia nearly 3000 years ago. In ancient Greece the physician Hippocrates had his patients throw them back and forth for injury prevention.

Reference	Description	Ø	Unit
06-030103	Mambo Max Medicine Ball - 1kg	20cm	1 pc
06-030104	Mambo Max Medicine Ball - 1,5kg	20cm	1 pc
06-030105	Mambo Max Medicine Ball - 2kg	20cm	1 pc
06-030107	Mambo Max Medicine Ball - 3kg	22cm	1 pc
06-030108	Mambo Max Medicine Ball - 4kg	22cm	1 pc
06-030109	Mambo Max Medicine Ball - 5kg	27cm	1 pc
06-000301	Mambo Max Medicine Ball rack (empty)		1 pc



01-RESISTIVE

02

111 Cb

06-STRE

Thumb Lock Wrist Weight <€

The Mambo Max Thumb Lock Wrist Weight enables you to fix the weight to your thumb, ensuring a fixed position during your entire workout. Being it boxing, running or Zumba, the weights will remain in place.



1111<u>ex</u>= (total of 300kg).

Available in: 2 - 4 - 6 - 8 - 10 - 12 - 16 - 20 - 24kg

120cm x 50cm x 70cm (l x w x h)

Reference	Weight	Description	Unit
06-050102	2kg	Mambo Max Kettlebell Sky Blue	1 pc
06-050104	4kg	Mambo Max Kettlebell Purple	1 pc
06-050106	6kg	Mambo Max Kettlebell Yellow	1 pc
06-050108	8kg	Mambo Max Kettlebell Red	1 pc
06-050110	10kg	Mambo Max Kettlebell Green	1 pc
06-050112	12kg	Mambo Max Kettlebell Blue	1 pc
06-050116	16kg	Mambo Max Kettlebell Black	1 pc
06-050120	20kg	Mambo Max Kettlebell Silver	1 pc
06-050124	24kg	Mambo Max Kettlebell Orange	1 pc
06-051001	Mambo N	Max Kettlebell Rack (empty)	1 pc
06-051101	Mambo N	Max Kettlebell Rack with 1 Kettlebell of each weight	set



28cm x 57cm x 183cm (Ixwxh)



The Mambo Max Kettlebells are an important part of any fitness or CrossFit training. Made from high quality iron surrounded with a colourful vinyl coating, they will not damage any surfaces when being stored. Starting from 2 kg up to 24 kg, we also offer you a complete set, including a premium quality Kettlebell Rack, which can hold up to 150kg per shelf



The Para Pro is a professional paraffin wax heating bath which is made to be used on 220/230 Volt power sockets and which can be used 24/7. In combination with our premium paraffin pearls it only needs between 2 and 3 hours to completely melt 2,5kg of paraffin. The temperature is adjustable from 40-60°C and is easy to control due to the digital thermometer.



- Paraffin Blocks:
- Green fully refined paraffin wax
- Melting point: 46°C
- Tea Tree scented




When ordering the Para Pro, you receive 2,5kg of white paraffin pearls or 6 blocks (2,72kg) of Tea Tree Green Paraffin free of charge.

Para Pro Specifications:

- Outer dimensions: 39cm x 25cm x 18cm (l x w x h)
- Inner dimensions: 32,5cm x 18cm x 13cm (l x w x h)
- Weight: 3,1kg
- 220 V: No transformer needed, which will lengthen the durability.
- Low power consumption once the paraffin has melted
- 24-hours a day fully operational
- Fast melting between 2 3 hours!
- Adjustable temperature (40 60°C) ٠
- Digital thermometer

Paraffin Pearls:

- White fully refined paraffin wax
- Congealing point: 53°C
- Melting point: +/- 54°C
- Meets the purity requirements of BGA1, Ph.EUR.2 and FDA3
- 1. Bundesgesundheitsamt recommendation No. XXV B II; 2. European Pharmacopoeia, current edition; 3. Food and Drug Administration, § 172.886 Petrole



Reference	Description	Unit
07-100101	MoVeS Para Pro Bath + 2,5kg Paraffin Pearls - White EU-plug	1 pc
07-100101 UK	MoVeS Para Pro Bath + 2,5kg Paraffin Pearls - White UK-plug	1 pc
07-100206	MoVeS Para Pro Bath + Set of 6 Paraffin Blocks 46°C (2,72 kg) - Tea Tree Green EU-plug	1 pc
07-100206 UK	MoVeS Para Pro Bath + Set of 6 Paraffin Blocks 46°C (2,72 kg) - Tea Tree Green UK-plug	1 pc
07-110102	MoVeS Paraffin Pearls - 2,5kg - White	1 pc
07-110125	MoVeS Paraffin Pearls - 25kg - White	1 pc
07-110206	MoVeS Paraffin Blocks 46° Set of 6 x 453g (2,72 kg) - Tea Tree Green	1 pc
07-109801	MoVeS Insulating Gloves for Para Pro Bath Blue	1 pair
07-109802	MoVeS Disposable Gloves for Para Pro Bath	10 pairs

Comes with UK or EU plug!



PearlPack MD (€

PearlPacks are reusable Hot/Cold Packs filled with small pearls. Hot and Cold Therapy can be applied to relieve pain and inflammations associated with muscle strain and sprain. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The small pearls absorb heat or cold, which they can retain for 20 minutes. During these 20 minutes, these packs will relieve a part of the pain.

Not only do they look good, the products are very convenient in use. PearlPacks won't feel wet nor drip after freezing or heating them, avoiding a mess. All materials used are non-toxic and odourless.

Besides relieving pain from an injury or muscle strain, hot/cold packs are often used to reduce swelling after an operation. PearlPacks truly are multifunctional products.

PearlPack Knee/Elbow



This offers a solution for Hot/Cold therapy on the knee or elbow joints. The large opening ensures that your knee or elbow remains free and can continue to move freely, while the Pearl technology sends direct relief to the pain. The flexible velcro strap keeps the pack firmly in place.

PearlPack Eye Mask

swollen eyes or sinusitis.

PearlPack Back Wrap



07-041101

07-040201



This large pack divided in 3 segments is equipped with velcro strap. It fits perfectly around the patient's waist for an optimal pain relief on the back. The pearls can be pushed from one segment to the other if one area of the back needs less or more thermal treatment.

PearlPack Sports Pack + strap



The Pearl Pack Sport pack is updated and now comes with a velcro strap.



Reference	Description	
07-040102	MoVeS PearlPack - Sport Pack + strap	UPDATED
07-040201	MoVeS PearlPack - Back Wrap	
07-041101	MoVeS PearlPack - Eye Mask	
07-040301	MoVes PearlPack - Knee/Elbow	

01-RESISTIVE

02-HAND

03-UPPER

20

LOWER BODY

04-WELLNESS

05-FITNESS

07-THERMAL

58



01-RESISTIVE

02-HAND

BODY

03-UPPER & LOWER

04-WELLNESS

05-FITNESS

06-STRENGTH

07-THERMAL



This mask adapts to the patient's face and has a fabric strap. Hot treatment can be applied when eyes feel dry or if there's an infection of the eyelid. Cold treatment will help when suffering from headaches,




















MoVeS Roll-On packs combine hot/cold treatment with practicality. The tubular shape and flexible material allow an easy roll-on application, even after being heated or frozen. The various sizes will fit perfectly over finger, wrist, ankle or knee. This design ensures that the whole area will benefit from the hot or cold treatment.

From the same factory as the original Torex[®] products. Reusable, non-toxic and latex free. Available in 5 sizes.





Reference	Description	Inside Ø	Circumference	Uses:	Unit
07-050101	MoVeS Roll-On Pack Finger	2,5cm	One size fits all	Finger	1 pc
07-050102	MoVeS Roll-On Pack Small	10,0cm	Fits up to 25cm	Hand, Wrist, Arm, Elbow	1 pc
07-050103	MoVeS Roll-On Pack Medium	12,5cm	Fits 25cm - 38cm	Hand, Wrist, Arm, Elbow, Ankle	1 pc
07-050104	MoVeS Roll-On Pack Large	15,0cm	Fits 38cm - 53cm	Elbow, Ankle, Foot, Calf, Knee	1 pc
07-050105	MoVeS Roll-On Pack Extra Large	17,5cm	Fits 53cm - 58cm	Knee, Thigh	1 pc

Reference Description

07-030201 Torex Hot/Cold Pro Ankle Boot 1 pc

Torex Pro Ankle Boot \in

Torex[®] Pro Ankle Boot A high performance ankle hot/cold and compression brace designed by PRO Athletic Trainers.

Unit





Hot/Cold Pack MD C€

MoVeS The MoVeS Hot/Cold Packs are usable for hot as well as cold therapy. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The versatile, reusable packs deliver hot and cold therapy effectively and economically. The use of the MoVeS Hot/Cold Pack is unlimited and its material is non-toxic. Even below -20°C it remains soft. The MoVeS Hot/Cold Pack can be bent in any angle, allowing maximum coverage of joint and muscle. All Hot/Cold Packs can be heated in the microwave and can be frozen in the freezer. Temperature holds for approximately 30 minutes.

MoVeS Hot/Cold Pack | Classic - Gel

MoVeS Hot/Cold Pack | Standard - Gel MoVeS Hot/Cold Pack | Soft Touch - Clay

The classic (and economical) transparent Hot/Cold Pack. The special gel is non-toxic and biodegradable.

Our standard Hot/Cold Pack is a nylon covered gelpack. The special gel is non-toxic and biodegradable.

Sold per 10 pieces

MoVeS

OT / COI

PACK



Protection sleeve

or home use only

100% polyester Protection Sleeve for hot/cold packs. Machine washable. Available in 3 sizes.

07-010206: XXL Hot/Cold Pack

This is ideal for pain relief on shoulders and back and has a nylon cover.



Moves HOT / COLD PACK

Dimension			MoVeS Hot/Cold	
		Classic (gel)	Standard (ge	
	15cm x 25cm	Small	07-010102	07-010202
	20cm x 30cm	Medium	07-010103	07-010203
	25cm x 35cm	Large	07-010104	07-010204
	33cm x 47cm	XXL	n/a	07-010206
	40cm x 20cm	Tri-Sectional	n/a	07-010211

02-HAND

01-RESISTIVE

05.

02-HAND

03-UPPER & LOWER BODY

04-WELLNESS



The Soft Touch is improved with a more user friendly surface. The frost-free surface allows usage in direct contact with the skin. Filled with clay, a 100% natural product.

Sold per 10 pieces







Ear & Forehead Thermometer

The infrared thermometer is used for measuring body temperature by two-in-one functions: Ear and Forehead

- temperature. It is suitable for medical unit and home use. Fast and easy temperature measurement by advanced infrared technology
- Range between 34.0°C~43.0°C ٠ Accuracy; $<35.0^{\circ}C$ and $>42.0^{\circ}C = \pm 0.3^{\circ}C$ 35.0°C~42.0°C = ±0.2°C
- Tri-colour backlit;
 - <37.4°C = Green
 - 37.4°C~38.0°C = Yellow
 - >38.0°C = Red
- Memory mode ٠
- Suitable for both infants and adults ٠
- Easy operating and cleaning
- Automatic power off after 60 seconds
- 2 AAA alkaline batteries





Reference	Description	Unit
08-070301	Infrared Thermometer Ear & Forehead auto-sensing 2 x AAA	1 pc
08-070303	Infrared Thermometer Shotgun Forehead auto-sensing 2 x AAA	1 pc










Blood Pressure Monitor MD CE

The blood pressure monitor is used for measuring blood pressure and pulse rate. It is suitable for people older than 12 years.



Reference	Description	Unit
08-070501	Blood Pressure Monitor + Cuff	1 pc

Fingertip Pulse Oximeter MD C€ cb m

A Pulse Oximeter is a very important and common device to check patient blood-oxygen saturation (SpO₂) level and pulse rate. As a small, compact, simple, reliable and durable physiological monitoring device, Fingertip Pulse Oximeter greatly enhances patient care.





•

08-070201

٠

Junior Oximeter

Rechargeable oximeter with OLED-screen, specifically designed for kids from 3-12 years old.

- Rechargeable via USB-cable •
- Dual colour OLED displays SpO₂, PR, pulse bar and waveform
- 4-direction & 6-mode display for convenient readings •
 - Set-up alarm range of SpO, and pulse rate
- Menu-function setting (beep sounds, etc.)
- Lithium battery; automatically power off ٠
- Gravity function, automatic rotating reading option

Reference	Description
08-070201	Fingertip Pulse Oximeter Adult OLED 2 x AAA
08-070202	Fingertip Pulse Oximeter Junior 3-12y old OLED

02-HAND

٠

.

62

07-THERMAL





Adult Oximeter

Features:

- Dual colour OLED displays SpO₂, PR, pulse bar and waveform
- SpO₂; measurement range between 70%~99%
- Accuracy; ±2% on the stage of 80%~99%
- PR; measurement range between 30BPM~240BPM
- Accuracy; ±1BPM or ±1% (the larger one)
- Low power indicator
- 6 display modes with adjustable brightness
- Automatically power off
- AAA alkaline batteries

	Unit
	1 pc
echargeable	1 pc

08-070202

02-HAND

08-EVALUATION

W Hydraulic 3-piece hand evaluation kit \in

SAEHAN

Portable carrying case contains all instruments necessary to evaluate the strength and range-of-motion of the hand and fingers.

Set includes:

- hydraulic hand dynamometer •
- hydraulic pinch gauge
- stainless steel finger goniometer 15cm 180° per 5°



01-RESISTIVE

02

HAND

03

Reference	Description	Unit
08-010113	Saehan Hydraulic Hand Evaluation Kit 3 pieces	Set



Portable carrying case contains all instruments necessary to evaluate the strength and range-of-motion of the hand and fingers.

SAEHAN

Set includes:

- hydraulic hand dynamometer
- hydraulic pinch gauge
- stainless steel finger goniometer 15cm 180° per 5°
- 2-point discriminator with 3rd point
- Wartenberg pinwheel
- finger gauge
- functional finger motion gauge



Unit

Set

Hydraulic hand dynamometer <€



hard carrying case.

Range: 0 - 90kg (200lb)

Reference	Description	Uni
08-010101	Saehan Hydraulic Hand Dynamometer	1 p

Digital hand dynamometer with Data Transfer Software- DHD-1 $\subset \in$



Designed for reliable grip-force evaluations from 0 to 90kg (200lb). It has an easy-to-read LCD display. This hand dynamometer features an electronic zero calibration system and up to 40 test data can be stored in memory. It is powered by a rechargeable Li-Ion polymer battery and data will not be lost when changing the battery. The Data Transfer

transfer the test results into .txt-files, which can be easily accessed using Notepad or MS-Word for editing and printing.

This Digital Hand Dynamometer hard carrying case includes a Digital Hand Dynamometer, the DTS-software and a USB cable.

Range: 0 - 90kg (200lb)

Reference Description

08-010201 Saehan Digital Hand Dynamometer - DHD-1

Digital hand dynamometer with G-STAR[™] Software - DHD-3

The DHD-3 Digital Hand Dynamometer (SH1003) prosame features as the DHD-1 when used vides the a stand-alone device, but adds automated as Rapid Exchange Grip Testing. Using Auto Mode, up to 10 strength tests for each hand as well as the time allowed to switch hands can be configured to automate the testing procedure.



The Grip Strength Testing and Research (G-STAR) is an easy to use software tool which automates grip strength testing and endurance testing in both clinical and research environments.

Range: 0 - 90kg (200lb) / Minimum requirements: Windows XP + USB port

Reference Description 08-010202 Saehan Digital Hand Dynamometer - DHD-3 (incl.

Reference Description 08-010117 Saehan Hydraulic Hand Evaluation Kit | 7 pieces

64

08-EVALUATION

07

THERMAL


SAEHAN

The Hydraulic Hand Dynamometer gives accurate grip strength readings without the patient being able to "feel" the handle move. The adjustable handle has 5 positions to accommodate to any hand size. The indicator remains at the patient's maximum reading until it is reset. Has a scale in kilograms. Comes in a





Unit
1 pc



	Unit
. G-STAR software)	1 pc

Smedley Hand Dynamometer <€



Push-Pull Dynamometer MD C€

A simple, easy-to-use, ergonomically designed instrument that objectively measures push, pull and lift forces for manual muscle testing, functional capacity evaluation and job task evaluation at a remarkably affordable price. Because the instrument is lightweight, small and portable, you can perform precise, objective evaluations in your practice or at the patient's location. The easy-toattach dual grip handle can be used when measuring larger forces.

Lift platform and chain to perform lifting evaluation not included.

\rightarrow	Dual Grip
	Range: 0 - 45kg (100lb)

Reference	Description	Unit
08-010501	Baseline Hydraulic Push-Pull Dynamometer	1 pc
08-010502	Dual Grip for Push-Pull Dynamometer	1 pc

Body Fat Scale MD C€

Baseline body fat scale is made with tempered glass. This scale measures weight, body fat and water % and memorises data up to 8 persons. Measurements range from 4% - 45% body fat; 37.8% - 66% body water and 330lb/150kg weight capacity. The required 1 x 3V battery is included.

Resolution: - 0.1% bodyfat - 0.1% water

Dimensions: 35cm x 35cm x 6cm Product weight: 2kg

Reference	Description	Unit
08-070101	Baseline BMI Body Fat Scale	1 pc





SAEHAN



08-020102: Mechanical Pinch Gauge

Used for hand-strength testing by measuring tip, key and palmar pinch force. Made of durable painted stainless steel. Mechanical properties will not degrade over time. It can serve as a diagnostic tool or for monitoring the progress of rehabilitation.



Reference	Description	Unit
08-020101	Saehan Hydraulic Pinch Gauge	1 pc
08-020102	Saehan Mechanical Pinch Gauge	1 pc

Squeeze Dynamometer \in



Reference	Description	Uni
08-010301	Saehan Squeeze Dynamometer 70kg	1 pc

Back-Leg-Chest Dynamometer <€

Measure strength of back, leg and chest. Base provides secure footing. Chain length is adjustable to accommodate for height differences or to vary the point of force application. Shows pounds and kilograms. Pointer remains at maximum until reset. Comes fully assembled with 1,5m chain. Range: 0 - 300kg (660lb)

Reference	Description	Unit
08-010601	Saehan Back-Leg-Chest Dynamometer	1 pc

Pinch Gauge \in

02

04-WELLNESS

07-THERMAL

08-EVALUATION





08-020101: Hydraulic Pinch Gauge

Uses the improved hydraulic system of the hand dynamometer to assure convenience, product reliability, measurement accuracy and repeatability. This results in a more accurate test result for all pinch tests.









The Squeeze (bulb) Dynamometer is the inexpensive device to reliably measure hand and finger strength. The dynamometer has an indicator dial that remains at the maximum reading until reset. Unit is calibrated in pounds and kilograms.








Fat Caliper <€



08-040101

Professional light-weight aluminium caliper. Measurement arms with floating tips to assure accurate and reproducible results. 65mm scale is large and easy to read. Corrective scale for sex and age included in multilingual manual.

- Professional, very accurate unit with spring loaded levers and low-pressure bearings for constant pressure.
- Scale permits reading up to 65mm
- Aluminium casted to insure years of reliable use.
- Includes carrying case ٠





Measurement Tape with Hands-free Attachment

SAEHAN

The Plastic Skinfold Caliper from Baseline is the inexpensive alternative for skinfold analysis. Spring loaded arm assures reproducible "pinch". Read directly from the 50mm scale. Ideal instrument for screening large population. The Plastic Skinfold Caliper comes with a booklet and tables.



The FatTrack® FITNESS TK Body Transformation Toolkit enables you to measure (and track) Body Fat %, Lean Body Mass, your workouts, macronutrient information and more by yourself - Easily, accurately and in the privacy of your own home!

Each FatTrack FITNESS TK Includes:

- Updated FatTrack[®] Fitness Digital Body Fat Caliper
- User Manual with instructions
- Body Tracker Software
- MyoTape Body Tape Measure
- FREE copy of "Improve Your Measuring Skill" with practical pointers for getting the best measurement readings and valuable information on ideal weight and muscle mass.

Reference	Description	Unit
08-040101	Saehan Skinfold Fat Caliper Aluminium - Green	1 pc
08-040104	Baseline Skinfold Fat Caliper Plastic - Black	1 pc
08-040102	FatTrack Skinfold Fat Caliper Digital	1 pc

Measurement Tape MD C€



Reference Description Unit Baseline Measurement Tape with Gulick Attachment - 150cm 08-080101 1 pc 08-080102 Baseline Measurement Tape with Gulick Attachment - 180cm 1 pc 08-080103 Baseline Measurement Tape with Hands-free Attachment - 150cm 1 pc 08-080201 Baseline Finger Circumference Gauge - 15cm 1 pc

 \sim

Tuning Fork

These aluminium Tuning Forks are developed to measure sensitivity to vibration which can be an indicator of nerve damage.

It is also used for foot exam by diabetic people to check for signs of neuropathy.

	Reference	Description	Unit
_	08-050301	Tuning Fork 128cps	1 pc
	08-050302	Tuning Fork 256cps	1 pc

Inclinometer

Baseline's Bubble and Digital Inclinometers are used to measure range-of-motion (ROM).

Some neck and back measurement protocols require the simultaneous use of 2 inclinometers.



Reference	Description	Unit
08-060101	Baseline Bubble Inclinometer	1 pc
08-060201	Baseline Digital Inclinometer	1 pc

Scoliometer

Scoliometers are inclinometers designed to measure trunk asymmetry or Angle of Truncal Rotation (ATR) used in scoliosis screening programmes. They are small, non-invasive devices which are placed over the spine while a person is in a forward bending position. The scoliometer reading indicates degree of trunk asymmetry (rotation). Research shows that most people have a small degree of spinal rotation. However, depending on the measurement, spinal rotation is a possible indicator that scoliosis may be present.

Baseline Plastic Scoliometer



Reference	Description
08-060102	Pedi-Scoliometer
08-060104	Baseline Plastic Scoliometer



01-RESISTIVE

02-

-HAND

03-UPPER & LOWER BODY

04-WELLNESS

68





MD (E



MD (E

Pedi-Scoliometer



Unit	
1 pc	
1 pc	

01-RESISTIVE

02-HAND

03-UPPER & LOWER

04-WELLNESS

05-FITNESS

06-STRENGTH

07-THERMAL



Reference	Description	Unit
08-030401	Baseline Plastic Absolute+Axis Goniometer - 30cm 360°	1 pc
08-030402	Baseline Plastic Absolute+Axis Goniometer - HiRes™ - 30cm 360°	1 pc
08-030411	Baseline Plastic Absolute+Axis - Attachment Only for 30cm Goniometer	1 pc
08-030601	Baseline Plastic Finger Goniometer - 9cm Flexion to Hyper-Extension	1 pc

Lafayette (€



Reference	Description	Unit
08-030701	Lafayette Gollehon Extendable Goniometer with Magnifying Lens - 24cm - 66cm 180°	1 pc

Goniometer Kit <€



The Plastic set includes:

- Plastic Goniometer 15 cm | 0° to 360° per 1°
- Plastic Goniometer 20 cm | 0° to 360° per 1°
- Plastic Goniometer 30 cm | 0° to 360° per 1° ٠
- Rulong Goniometer 18 cm | 0° to 360° per 2°
- Pocket Goniometer 15 cm | 180° x2
- Plastic Flexion-Hyperextension goniometer

Ken Sten St

SAEHAN

Quality Soft Pouch Case

The Stainless Steel set includes:

- Robinson pocket goniometer 15cm 180°
- Stainless steel goniometer 18cm 180°
- Stainless steel goniometer 35cm 180° x2
- Stainless steel goniometer 35cm 360°
- Stainless steel finger goniometer 15cm
- X-ray goniometer 20cm 180°





)° per 5°

30° per 5°

1 pc

1 pc



Reference	Description
08-030301	Saehan Finger Goniometer 9cm Stainless Steel 180
08-030302	Saehan Finger Goniometer 15cm Stainless Steel 18

Reference	Description	Unit
08-030211	Saehan Goniometer Set Stainless Steel 6 pieces	Set
08-030121	Saehan Goniometer Set Plastic 6 pieces	Set

02

HAND

0

04-WELLNESS

08-EVALUATION

70

71

Monofilament MD (E



Reference

08-050401

08-050203

Touch-Test 2-point Discriminator


08-050206



08-050401

Description

_	Neurological	Reflex	Instrum



Reference	Description	Unit
08-050203	Touch-Test 2-point Discriminator	1 pc
08-050204	Baseline 2-point Discriminator	1 pc
08-050205	Baseline 3-point Discriminator	1 pc
08-050206	Baseline 2-point Discrim-A-Gon, 2 Disk Set	1 pc

05

-FITNESS

01-RESISTIVE

02

-HAND

03-

UPPER & LOWER BODY

04-WELLNESS

07-THERMAL

72









in a diagnostic evaluation. It features a double-sided head composed of two rounded rubber surfaces for

you sit with your feet against the "sit and reach box" and flex forward moving the slide bar. Flexion is read directly from the calibrated bar.



01-RESISTIVE

Cb m

01-RESISTIVE

02-HAND

03-UPPER & LOWER BODY

04-WELLNESS

TENS device MD (E



The Neuropex TENS is accepted to be far more accurate than a standard TENS machine.

- Highly user friendly with 3 programs and modes:
 - B (burst)
 - C (continuous)
 - M (Pulse Rate Modulation)
- Regulated pulse width control for adjustable intensity and stimulation strength
- Clear plastic window allows controls to be viewed during ٠ operation
- Easy to use and very compact ٠
- Adjustable intensity and stimulation strength
- Protective case for easy electrode storage and transportation •
- Complete set of electrodes, in-built timer, complimentary spare • battery







Neuropex Long Lasting Self-Adhesive Electrodes ٠

Resealable pouch for safe storage 4 Electrodes per pack

Interferential devices

- 50mm x 50mm
- Ideal for use with all Neuropex TENS and



О

0

We have nice products, but you would like your own brand or your customer's name on the band, loop, superloop ...? We'll do it! Together with our manufacturers from all the products in our range, we can create your desired product with the look you want. We can replace our logo with yours, change the colour, even the packaging can be tweaked to your needs. Private Labels will allow you to design the product the way you want to match it to your company's or customer's identity.



Private Label products have the same quality as all our products, but are putting your brand in the spotlights. This is perfect if you need a homogeneous range of products or want to increase your brand awareness. Or maybe you're looking for a one time opportunity: a promotional product during your company's special occasion, a gift for visitors or an event you're supporting?

Please contact us for more information about pricing and volume.

Product descri



Unit Reference Description 08-070401 2-Channel Professional TENS unit 1 pc 08-070411 Replacement TENS Electrode - 4 piece set set

74

PRIVATE LABEL



Product description	Minimum Order Quantity
Resistive Exercise Band	5.000m per colour
Resistive Exercise Loop	500 pcs per colour
Soft-Over-Ball	3.000 pcs per colour
Hollow Foam Roller	300 pcs
AB Gym Ball	500 pcs
NBR Comfort Gym Mat NBR Xtra Comfort Gym Mat	300 pcs per colour
Pearl Packs	500 pcs
Plastic Goniometer	2.000 pcs
Various other products	Contact us for more info

1000

Ó

AFTER-SALES ///



Measurement Instruments are calibrated at manufacturing with a full one year guarantee, before sending them to users all over the world. Ideally, these instruments should be recalibrated every year.

In order to make this happen in an easy and inexpensive way, Hydraulic Hand Dynamometers (Saehan, Jamar, Baseline...) and Hydraulic Pinch Gauges are recalibrated by our experts at the MVS In Motion Centre in Tisselt - Belgium. This operation consists of re-establishing the exact relationship between the force applied to a dynamometer and the measurement it shows. Small repairs can be done as well.



Reference	Description	Unit
08-019901	Calibration Unit	1 pc
08-019900	Calibration & Repair	-

MVS IN MOTION REPAIR CENTRE



Since machines are designed and assembled by humans, they are not flawless. That is why MVS In Motion installed an inhouse repair centre, to swiftly and correctly solve any problems that may occur with our products.

No matter what the cause, be it faulty fabrication or misuse, our technician will shine his light on the problem.

Send an email to our customer service department to inform us on the type of problem and the serial number/ production date of the device. After our authorisation, please send over the patient, accompanied by your name and client number and a brief description of the symptoms to info@mvs-in-motion.com. We'll do our best to have him cured in no time.

3-in-1 Foam Roller	43	Flex-Ion	25	Ring Tubes	20
Accessories	24	Foam Hand Grip	32	Roll-On Packs	60
Adjustable Hand Grip	32	Foot Pump	48	Saehan Cycle	34
Aerobic Step	48	Goniometers	70-71	Safety Handles	24
Ankleciser	21	Gym Ball	51	Scoliometer	69
Babinski Hammer	73	Gym Mat	47	Shoulder Rope Pulleys	22,23
Back Wrap	59	Hand Pump	48	Shoulder Tube Pulley	22
Back-Leg-Chest Dynamometer	67	Handmaster Plus	27	Sit & Reach Unit	73
Balance Boards	50,51	Hollow Foam Roller	43	Sitting Cushions	44
Balance Pads	53	Hot/Cold Packs	61	Smedley Hand Dynamometer	66
Balance Trainers	51,52	Hydraulic Hand Dynamometer	65	Soft Handles	24
Ball Accessories	48	Hydraulic Hand Evaluation Kits	64	SoftMeds	55
Ball Stabiliser	48	Incline Boards	39	Soft-Over-Ball	45
Ball Stacker	48	Inclinometer	69	Spiky Foam Roller	43
Band	8-11,17	Kettlebells	57	Sport Pack + Strap	59
Band & Tubing Scissors	24	Kinesiology Tape	41	Squeeze Balls & Eggs	31
Band Latex-Free	13	Knee Rope Pulley	22	Squeeze Dynamometer	67
Bar	20	Knee/Elbow Wrap	59	Stepping Stones	42
Blood Pressure Monitor	63	Loop	15,19	Superloop	19
Body Fat Scale	66	Massage Ball, Roll & Ring	46	Swim Noodle	54
Body Twister	52	Massage Stick	43	Taylor Hammer	73
Buck Hammer	73	Massage Tools	46	TENS Device	74
Comfort Putty	33	Measurement Tape	68	Theraflex Putty	33
Constant Force X-Tend	27	Medicine Ball	57	Thermometers	62
Cuff-Ring Tube Loop	20	Monofilament	72	Thumb Lock Wrist Weight	56
Curamotion	38	Multi-Grip Exerciser	21	Torex Pro Ankle Boot	60
Dejerine Hammer	73	Needle Pump	48	Triggerpoint Rollers	43
Digi-Extend	26	Neurological Instruments	73	Troemner Hammer	73
Digi-Flex	26	O-Ring Tube Loop	20	Tube	14,18
Digital Hand Dynamometers	65	Oximeters	63	Universal Assist	24
Discriminators	72	Oxycycles	39-37	Varigrip Sport	28
Door Anchor	23	Paraffin Bath	58	Varigrip Therapy	28
Dumbbells	55	Peanut Ball	51	Varigrip Uno	28
Dynadome	49	PearlPacks	59	Wall Mount	24
Ekilibres	52	Pedal Exercisers	39	Wartenberg Pinwheel	73
Electric Compressor Pump	48	Pilates Foam Roller	45	Wheeleo	35
Ergonomic Foam Roller	43	Pilates Ring	45	Wide Loop	15
Eye Mask	59	Pinch Gauge	67	Work-Out Station	52
Fabric Elastic Band	49	Power-Webs	29-30	Wrist & Ankle Weight	56
Faster Blaster Pump	48	Prohands Series	26	Wrist Exerciser	31
Fat Caliper	68	Pumps	48	Yoga Block	42
Figure 8-Ring Tube Loop	20	Push-Pull Dynamometer	66	Yoga Mat	42
Fit Stretch	34	Putties	33	Yoga Strap	42

<u>INDEX</u>

///

TERMS & CONDITIONS

- 1. Terms: Prepayment, unless otherwise specified. Buyer should bear all bank charges.
- 2. Minimum purchase value: An administration fee of €10,00 will be added on all orders below the minimum purchase value of €250,00.
- 3. Freight Terms: Prices are ex-works Tisselt, Belgium. All shipping/delivery costs are at charge of the buyer, unless otherwise specified. Visible damage to goods delivered, not duly noted on the consignment note of the carrier, cannot be compensated in any form.
- 4. Prices: All prices for merchandise purchased are based on seller's price list in effect on the time of the order. Prices can be subject to change without notice.
- 5. Claims: Unless otherwise provided for herein, claims of any kind in goods sold by seller shall be made by buyer to seller, in writing, within a period of seven days after such goods are delivered to buyer.
- 6. Credit and Exchange: Written authorisation must be obtained before returning merchandise for credit or exchange! Prior to return, buyer shall furnish seller an itemized inventory of the products to be returned, and any other information seller may request, such as product lot number, manufacturing date, expiration date, ... In no event shall seller accept products for return more than thirty (30) days after delivery by seller to buyer and products which are damaged, used or not in their original and unopened packaging. Seller will charge a handling fee of ten percent (10%) of the returned goods.
- 7. Finance Charge: A finance charge of 1% per month (12 % per annum) will be added to all amounts past due. A penalty of 15% on value, with a minimum of €40,00, for late payment will be added. Terms of payment are from date seller's invoice and not contingent upon delivery.
- 8. Governing Law: All disputes arising out of all sales by seller to buyer shall be governed by the laws of Belgium and place of jurisdiction is Antwerp, section Mechelen.
- 9. Warranties: Seller warrants that it has good title to the goods purchased by buyer and that its goods will conform to buyer's specifications, if any, except as otherwise provided herein. All other warranties, whether express or implied, and including, but not limited to, the warranty of merchantability and fitness for a particular purpose are excluded. In all events, buyer's sole and exclusive remedy shall be limited to the amount of the purchase price of any product on which a claim is made as its exclusive remedy for breach of the warranty expressed herein and for negligence and otherwise. Seller reserves the right to repair or replace warranted products in its descretion. Under no circumstances shall seller be liable to buyer or any other person for consequential or incidental damages.
- 10. Property: Seller has the exclusive property on the delivered products until buyer has paid the full amount due. This also implies that if buyer does not pay in time, seller has the right to, without prenotice, claim back the delivered goods.








The choice of physiotherapists for 50 years, with know-how and expertise at your service.

www.ferrox.it info@ferrox.it • +39 (0)4 38 77 70 99



VISIT OUR OTHER WEBSITES FOR

- PRODUCT INFORMATION
 EXERCISE CHARTS AND VIDEOS
- SHOP LOCATOR

WWW.MOVES-YOU.COM WWW.MOVES-BAND.COM WWW.MAMBOMAX.COM

MVS In Motion | Westdijk 150, 2830 Tisselt (Willebroek) | Belgium - Europe [T] +32 3 866 40 76 | www.mvs-in-motion.com | info@mvs-in-motion.com